

Suggested protocol for loading a DNA Ladder/marker

New England Biolabs

Abstract

This is the suggested protocol for use with λ DNA-Mono Cut Mix ([N3019](#)), ϕ X174 DNA-HaeIII Digest ([N3026](#)), pBR322 DNA-BstNI Digest ([N3031](#)), pBR322 DNA-MspI Digest ([N3032](#)), 2-Log DNA Ladder (0.1-10.0 kb) ([N3200](#)), 100 bp DNA Ladder ([N3231](#)), 1 kb DNA Ladder ([N3232](#)), Low Molecular Weight DNA Ladder ([N3233](#)), and 50 bp DNA Ladder ([N3236](#))

Citation: New England Biolabs Suggested protocol for loading a DNA Ladder/marker. **protocols.io**
dx.doi.org/10.17504/protocols.io.cq4vyv

Published: 31 Jan 2015

Before start

The following protocol is recommended for a 5 mm wide lane.

Protocol

Step 1.

Prepare loading mixture (**6 μ l** total volume):

✓ PROTOCOL

. [Loading Marker Mixture](#)

CONTACT: [New England Biolabs](#)

⊕ NOTES

New England Biolabs 31 Jan 2015

The components of the mixture should be scaled up or down, depending on the width of the agarose gel.

Step 1.1.

Distilled water, **4 μ l**

📄 AMOUNT

4 μ l Additional info:

Step 1.2.

6X Purple no-SDS Loading Dye, **1 μ l**

📄 AMOUNT

1 μ l Additional info:

🧴 REAGENTS

🐛 Gel Loading Dye, Purple (6X), no SDS - 4.0 ml [B7025S](#) by [New England Biolabs](#)

Step 1.3.

DNA Ladder, **1 μ l**

AMOUNT

1 µl Additional info:

Step 2.

Mix gently

Step 3.

Load onto the agarose gel

NOTES

New England Biolabs 31 Jan 2015

For long term storage, store at -20°C. If samples need to be diluted, use TE or other buffer of minimal ionic strength. DNA may denature if diluted in dH2O.