

# Cryopreserve cells

Bao Thai

## Abstract

Cryofreezing cells for preservation.

**Citation:** Bao Thai Cryopreserve cells. **protocols.io**

dx.doi.org/10.17504/protocols.io.mtsc6ne

**Published:** 25 Jan 2018

## Protocol

### Step 1.

Grow cells to 100% confluency in a 10cm dish, which can make 5 cryotubes with 1 mL each.

### Step 2.

Trypsinize cells (1.5 mL trypsin in a 10mm plate).

### Step 3.

Incubate at 37C for 5 minutes.

### Step 4.

Resuspend cells completely off the plate and add 5mL of DMEM media containing FBS. Transfer to a 14 mL tube.

### Step 5.

Centrifuge tubes at 100xg for 15 minutes.

### Step 6.

Decant supernatant carefully.

### Step 7.

Resuspend cell pellets completely in 5 mL of DMEM containing 5% DMSO.

### Step 8.

Add 1 mL of the solution containing cells into each cryogenic tube.

### Step 9.

Freeze cells slowly by putting the tubes in between tightly sealed styrofoam boxes at -80C for 2 days.

### Step 10.

Transferred tubes to liquid nitrogen for long term storage.