OPEN ACCESS



Quick Protocol for Monarch® Plasmid Miniprep Kit (NEB #T1010) Version 2

New England Biolabs

Abstract

This is the quick version of the Monarch® Plasmid DNA Miniprep Kit Protocol (NEB #T1010). For the full protocol, please click here.

Citation: New England Biolabs Quick Protocol for Monarch® Plasmid Miniprep Kit (NEB #T1010). protocols.io

dx.doi.org/10.17504/protocols.io.nytdfwn

Published: 21 Mar 2018

Guidelines

For detailed protocol and more information, visit www.neb.com/T1010

The full protocol is available <u>here</u>.

The video protocol is available here.

Before start

- All centrifugation steps should be carried out at 16,000 x g (~13,000 RPM).
- Add 4 volumes of ethanol (≥ 95%) to one volume of Plasmid Wash Buffer 2.
- If precipitate has formed in Lysis Buffer (B2), incubate at 30–37°C, inverting periodically to dissolve.
- Store Plasmid Neutralization Buffer (B3) at 4°C after opening, as it contains RNase A.

Materials

Monarch® Plasmid Miniprep Kit <u>T1010</u> by New England Biolabs

Protocol

Step 1.

Pellet 1–5 ml (not to exceed 15 OD units) bacterial culture by centrifugation for 30 seconds. Discard supernatant.

© DURATION

00:00:30

P NOTES

Danielle Freedman 21 Mar 2018

1.5 ml of culture is sufficient for most applications. Ensure cultures are not overgrown (12-16 hours is ideal).

Step 2.

Resuspend pellet in 200 µl Plasmid Resuspension Buffer (B1).

₽ NOTES

Danielle Freedman 21 Mar 2018

Vortex or pipet to ensure cells are completely resuspended. There should be no visible clumps.

Step 3.

Add 200 μ l Plasmid Lysis Buffer (B2), gently invert tube 5–6 times, and incubate at room temperature for 1 minute. Do not vortex.

O DURATION

00:01:00

NOTES

New England Biolabs 17 Feb 2016

Color should change to dark pink, and solution will become transparent and viscous. Do not vortex.

Step 4.

Add 400 μ l of Plasmid Neutralization Buffer (B3), gently invert tube until neutralized, and incubate at room temperature for 1 minute. Do not vortex.

O DURATION

00:02:00

P NOTES

New England Biolabs 17 Feb 2016

Sample is neutralized when color is uniformly yellow and precipitate forms. Do not vortex.

Step 5.

Centrifuge lysate for 2-5 minutes.

O DURATION

00:02:00

P NOTES

Danielle Freedman 21 Mar 2018

For culture volumes >1 ml, we recommend a 5 minute spin to ensure efficient RNA removal by RNase A. Pellet should be compact; spin longer if needed.

Step 6.

Carefully transfer supernatant to the spin column and centrifuge for 1 minute. Discard flow-through.

O DURATION

00:01:00

Step 7.

Re-insert column in the collection tube and add 200 μ l of Plasmid Wash Buffer 1. Centrifuge for 1 minute.

© DURATION

00:01:00

NOTES

Danielle Freedman 21 Mar 2018

Discarding the flow-through is optional.

Step 8.

Add 400 μ l of Plasmid Wash Buffer 2 and centrifuge for 1 minute.

O DURATION

00:01:00

Step 9.

Transfer column to a clean 1.5 ml microfuge tube.

NOTES

New England Biolabs 17 Feb 2016

Use care to ensure that the tip of the column does not come into contact with the flow-through. If there is any doubt, re-spin the column for 1 minute.

Step 10.

Add \geq 30 μ l DNA Elution Buffer to the center of the matrix. Wait for 1 minute, then spin for 1 minute to elute DNA.

NOTES

Danielle Freedman 21 Mar 2018

Nuclease-free water (pH 7–8.5) can also be used to elute the DNA. Yield may slightly increase if a larger volume of DNA Elution Buffer is used, but the DNA will be less concentrated. For larger size DNA, (≥ 10 kb), heating the elution buffer to 50°C prior to use can improve yield.