

# Non-Centrifuged Plaque Assay Sample Steps

Matthew Sullivan

## Abstract

For [Transcriptomics During One-Step Growth Curves for Cellulophaga Phages protocol](#).

**Citation:** Matthew Sullivan Non-Centrifuged Plaque Assay Sample Steps. **protocols.io**

dx.doi.org/10.17504/protocols.io.der3d5

**Published:** 02 Feb 2016

## Protocol

### Step 1.

Pipet 10  $\mu$ l from the flask into 990  $\mu$ l of MSM in a 1.5 ml tube

#### NOTES

**VERVE Team** 15 Jul 2015

You are diluting your sample 100x:  $10^{-2}$ .

### Step 2.

Vortex briefly

### Step 3.

Do serial dilutions to make  $10^{-4}$  and  $10^{-5}$  dilutions

### Step 4.

Plate 100  $\mu$ l each from the  $10^{-4}$  and  $10^{-5}$  dilutions