

# Protocol D: Preparation of BSA for quenching pretreatment

David Booth

## Abstract

**Citation:** David Booth Protocol D: Preparation of BSA for quenching pretreatment. **protocols.io**

dx.doi.org/10.17504/protocols.io.h8wb9xe

**Published:** 22 Jun 2018

## Protocol

Make a 10% (w/v) solution of BSA

### Step 1.

- Dissolve 2 g of bovine serum albumin, fraction V (Sigma 10735078001) in 20 ml of 50 mM HEPES-KOH pH 7.5, 200 mM NaCl, and 20% (v/v) glycerol
- 2. Sterile filter BSA solution through a 0.22 µm filter.
- 3. Prepare a small aliquot of a 100 fold dilution of BSA in 8 M Urea and measure the optical absorbance at 280 nm with a UV-Vis spectrophotometer. Determine the concentration of BSA where 1 absorbance unit is a concentration of 0.667 mg/ml BSA in a 1 cm pathlength.
- 4. Dilute BSA to a final concentration of 50 mg/ml with 50 mM HEPES-KOH pH 7.5, 200 mM NaCl, and 20% (v/v) glycerol that has been sterily filtered through a 0.22 µm filter.
- 5. Make 1 ml aliquots and store at -20°C.

Sterile filter BSA solution through a 0.22 µm filter

### Step 2.

Quantify BSA solution by UV-Vis spectrophotometry

### Step 3.

- Prepare a small aliquot of a 100 fold dilution of BSA in 8 M Urea
- Measure the optical absorbance at 280 nm with a UV-Vis spectrophotometer
- Determine the concentration of BSA where 1 absorbance unit is a concentration of 0.667 mg/ml BSA in a 1 cm pathlength.

Dilute BSA to a final concentration of 50 mg/ml

### Step 4.

- Dilute with 50 mM HEPES-KOH pH 7.5, 200 mM NaCl, and 20% (v/v) glycerol that has been sterily filtered through a 0.22 µm filter

Store at -20°C

## Step 5.

---