

# Yeast DAPI Staining

Alan Cone

## Abstract

Stains the nucleus of budding or fission yeast with DAPI in order to see it with a fluorescence microscope.

**Citation:** Alan Cone Yeast DAPI Staining. [protocols.io](https://protocols.io)

[dx.doi.org/10.17504/protocols.io.eigbcbw](https://dx.doi.org/10.17504/protocols.io.eigbcbw)

**Published:** 10 Feb 2016

## Protocol

### Step 1.

Grow up yeast in liquid medium overnight.

#### NOTES

**Alan Cone** 10 Feb 2016

The OD may not matter too much here, but something in the 0.8 - 2 range is probably ideal.

### Step 2.

Add 333  $\mu$ L (or 1 volume) yeast culture to a 1.5 mL microcentrifuge tube.

#### AMOUNT

333  $\mu$ L Additional info:

### Step 3.

Add 666  $\mu$ L (or 2 volume) of 100% Ethanol to the 1.5 mL microcentrifuge tube.

#### AMOUNT

666  $\mu$ L Additional info:

#### REAGENTS

Ethyl alcohol, Pure 200 proof, for molecular biology [E7023](#) by [Sigma Aldrich](#)

### Step 4.

Let the yeast ethanol mixture sit at room temperature for 30-60 minutes.

#### DURATION

00:30:00

### Step 5.

Spin down yeast cells for 1 minute at 2500 RPM.

#### DURATION

00:01:00

### Step 6.

Pour out the supernatant and resuspend the pellet in 1mL of 1 x PBS, then centrifuge for 1 minute at 2500 RPM.

#### AMOUNT

1 ml Additional info:

#### REAGENTS

✓ 1X PBS (Phosphate-buffered saline ) by Contributed by users

### Step 7.

Pour out the supernatant and resuspend the pellet in 200 µL of a 1 x PBS / 1:2000 Dilution DAPI mixture.

📄 AMOUNT

200 µl Additional info:

📄 PROTOCOL

### . [PBS / DAPI 1:2000 Dilution Mixture](#)

CONTACT: [Alan Cone](#)

#### Step 7.1.

Add 1 mL 1 x PBS to a 1.5 mL microcentrifuge tube.

📄 AMOUNT

1 ml Additional info:

📄 REAGENTS

✓ 1X PBS (Phosphate-buffered saline ) by Contributed by users

#### Step 7.2.

Add 0.5 µL of a 2.5 mg/mL (or a 1:2000 dilution) of DAPI to the 1 x PBS.

📄 REAGENTS

✓ DAPI (2.5mg/mL) by Contributed by users

### Step 8.

Add one drop of the yeast suspended in the PBS / DAPI mixture onto a microscope slide, add a coverslip on top, and go observe the stained yeast.

📌 NOTES

**Alan Cone** 10 Feb 2016

Make sure after you do this part to go look at it within a few hours of resuspending in PBS / DAPI. The sooner you can get to a microscope the better.

## Warnings

DAPI is light sensitive.

You must centrifuge at low RPM to avoid displacement of the nucleus.