



Cas9 RNP nucleofection for CD34+ HSPCs using Lonza 4D Nucleofector
Forked from [Cas9 RNP nucleofection for cell lines using Lonza 4D Nucleofector](#)

[Mark DeWitt](#)¹
¹IGI/UC Berkeley
[dx.doi.org/10.17504/protocols.io.hdnb25e](https://doi.org/10.17504/protocols.io.hdnb25e)

CornLab



Jacob Corn
ETH Zurich



PROTOCOL STATUS


Working
We use this protocol in our group and it is working

STEPS MATERIALS


NAME	CATALOG #	VENDOR
SFEM II		Stemcell Technologies
CC110		Stemcell Technologies
Mobilized Peripheral Blood CD34 Stem/Progenitor Cells	mPB015F	

Thaw HSPCs, plate at 300,000 cells/mL in stem cell expansion medium


- 1



SFEM II
by [Stemcell Technologies](#)



CC110
by [Stemcell Technologies](#)



Mobilized Peripheral Blood CD34
Stem/Progenitor Cells
Catalog #: [mPB015F](#)
- Count cells daily using hemacytometer
- 2Log counts for day 0 (thaw day), day 1, and day 2 (electroporation day).

Generally, the count will decline slightly from day 0 to day 1, and recover from day 1 to day 2. The total number of cells on day 2 will be slightly higher than on day 0.
- 3
- ✓ protocols.io
- 1
- 01/02/2019
- This is an open access protocol distributed under the terms of the Creative Commons Attribution License (<https://creativecommons.org/licenses/by/4.0/>), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Prepare RNP mix

4

5 Add Cas9 to sgRNA slowly while swirling pipette tip, should take 30s to 1 minute.

6

🕒 00:20:00

7 Count cells using hemacytometer.

8 For each nucleofection, pipette 200k cells into a 15 mL conical.

NOTE

For replicates, you can multiply the amount of cells and RNP mix as needed and mix in a single tube before electroporating separately.

Prepare Cells

9 Spin 200 x g for 10 minutes to pellet cells softly. While the cells are spinning, prepare plate and cuvette.

🕒 00:10:00

10 Place a 96 well plate in the tissue culture incubator.

Nucleofection

11 Prepare and label wells on 20uL nucleofection strips. Configure Lonza 4d to zap with code ER100.

12 Pipette off media from cells, gently but completely, using a P200. The pellet is very soft so be careful.

13 Resuspend cells in 20 μ L of P3 nucleofector solution using a P200.

14 Add the entire 10 μ L RNP mix to the 20 μ L resuspension and mix.

15 Add 1uL of 100uM donor DNA (100 pmoles) and mix well.

16 Add 27.5 μ L of cells/RNP electroporation mix to the multiwell cuvette, and cap. Pay attention to the orientation of the cap and cuvette in the nucleofector, which is noted in the manufacturer's instructions.

17 Insert cuvette into nucleofector and zap.

18 Immediately layer electroporated cells with 72.5 μ L of medium (SFEM II/Pen-strep/CC110) and rest for 5 minutes at room temperature.

19

Nucleofection

20 Pipette mixture out with a P200 into your pre-warmed 96 well plate. Wash each cuvette with an additional 100 μ L of medium, and add that as well (final plating volume 200 μ L).

21 Culture cells for several days in SFEM II/PenStrep/CC110 before genotyping by desired method. Consider including un-zapped controls to test viability.

Alternatively, transfer cells to expansion medium after 1-2 days. We have used StemCell Erythroid Expansion supplement for this purpose. Otherwise, the cells will not expand substantially if left in CC110



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited