

Orthopoxvirus real-time PCR

Judy Northill, David Warrilow, Ian Mackay

Abstract

A real-time PCR targeting the DNA-dependent RNA polymerase of Orthopoxviruses.

Citation: Judy Northill, David Warrilow, Ian Mackay Orthopoxvirus real-time PCR. **protocols.io**

dx.doi.org/10.17504/protocols.io.nxzdfp6

Published: 23 Mar 2018

Guidelines

- If using a different brand or model of real-time thermocycler, check the concentration of ROX is adequate.
- Method assumes the user is familiar with the thermocycler and software used to run the protocol.

Materials

SensiFAST™ Probe Lo-ROX Kit [BIO-84002](#) by [Bioline](#)

Protocol

Oligonucleotide sequences

Step 1.

| Name | 5'-3' |
|-------------|--|
| OPV2018-F | CGTACMGGAACTRGCTAGA |
| OPV2018-R | AGCGTATTACCTATAACTACTTGTCCGTA |
| OPV2018-FAM | 6FAM- ATCATTAAAAAGATGGAGGATATGGTGGTHGA -BHQ1 |

Step 2.



REAGENTS

SensiFAST™ Probe Lo-ROX Kit [BIO-84002](#) by [Bioline](#)

REACTION SET-UP

Step 3.

Assay has been used on both a Rotor-Gene 6000 / Rotor-Gene Q 5-plex using 100-place rotor discs.

Total reaction volume is 20µL.

Prepare sufficient for number of reaction plus a 'dead volume' usually 2 extra. Adjust as necessary if using a robotic dispenser.

| Reagent | Volume (uL) | Final concentration |
|----------------------------|-------------|---------------------|
| Nuclease-free water | 4.87 | |
| SensiFast Probe Lo-Rox mix | 10 | 1X |
| OPV2018-F | 0.05 | 500nM |
| OPV2018-R | 0.05 | 500nM |
| OPV2018-FAM | 0.03 | 150nM |
| TOTAL VOLUME | 15 | |

AMPLIFICATION

Step 4.

The assay has been optimised and validated for the Rotor-Gene 6000 and Rotor-Gene Q thermocyclers.

PCR

| | | |
|------|------|-----|
| 50°C | 5min | |
| 95°C | 2min | |
| 95°C | 3s | 40X |
| 60°C | 30s* | |

*Data acquisition step

RESULT ANALYSIS

Step 5.

The definition used for a satisfactory positive result from a real-time fluorogenic PCR should include each of the following:

1. A **sigmoidal curve** – the trace travels horizontally, curves upward, continues in an exponential rise and followed by a curve towards a horizontal plateau phase
2. A **suitable level of fluorescence** intensity as measured in comparison to a positive control (y-axis)
3. A defined threshold (C_T) value which the fluorescent curve has clearly exceeded (Fig.1 arrow) and which sits early in the log-linear phase and is <40 cycles
4. A flat or non-sigmoidal curve or a curve that crosses the threshold with a C_T value >40 cycles is considered a negative result

5. NTCs should not produce a curve




Figure 1. Examples of satisfactory sigmoidal amplification curve shape when considering an assay's fluorescent signal output. The crossing point or threshold cycle (C_T) is indicated (yellow arrow); it is the value at which fluorescence levels surpass a predefined (usually set during validation, or arbitrary) threshold level as shown in this normalized linear scale depiction. LP-log-linear phase of signal generated during the exponential part of the PCR amplification; TP-a slowing of the amplification and accompanying fluorescence signal marks the transition phase; PP-the plateau phase is reached when there is little or no increase in fluorescent signal despite continued cycling.