

# QA-Bio PNGase F Non-denatured Glycoprotein Protocol

Mike Gibson

## Abstract

PNGase F cleaves N-linked (asparagine-linked) oligosaccharides from glycoproteins. The enzyme deaminates asparagine to aspartic acid, leaving the oligosaccharides intact. Denaturation increases the rate of cleavage. Most native proteins can still be completely N-deglycosylated but incubation time must be increased. The enzyme will remain fully active under reaction conditions (37°C) for at least 96 hours. PNGase F will not remove oligosaccharides containing Alpha-(1,3)-linked core fucose commonly found on plant glycoproteins; for this purpose, use peptide N-glycosidase A.

**Citation:** Mike Gibson QA-Bio PNGase F Non-denatured Glycoprotein Protocol. **protocols.io**  
dx.doi.org/10.17504/protocols.io.seaebae

**Published:** 06 Aug 2018

## Materials

 PNGase F [E-PNG01](#) by [QA-Bio Inc](#)

## Protocol

### Step 1.

Add up to 200 µg of glycoprotein to an Eppendorf tube.

### Step 2.

Adjust to 40 µl final volume with de-ionized water.

### Step 3.

Add 10 µl 5x Reaction Buffer 7.5

### Step 4.

Add 2.0 µl of enzyme to the reaction.

### Step 5.

Incubate 3-24 hours at 37°C.