

untitled protocol

Yuan Yao

Abstract

Citation: Yuan Yao untitled protocol. **protocols.io**

dx.doi.org/10.17504/protocols.io.qaddsa6

Published: 22 May 2018

Protocol

Step 1.

Incubate E.coli cells 1:100 in 10 ml LB medium with antibiotics.

Step 2.

OD600=0.4 , collect all cells at 4 °C, 8000 rpm for 10 min.

Step 3.

Wash with 1 ml 1*PBS one time.

Step 4.

Resuspend in 500 ul 1 * PBS and sonicate.

Step 5.

Spin at 10000 rpm for 15 min.

Step 6.

Transfer supernatant into a new 1.5 ml tube and sit on ice.

Step 7.

Load into 120% SDS-PAGE gel , run at 200v for 45 min.

Step 8.

Get the PVDF membrane and filter paper
PVDF: Put in methanol (5s) then in transfer buffer
Filter: transfer buffer
Gel: transfer buffer

Step 9.

Sit on Room Temperature for 20 min.

Step 10.

Transfer membrane at 25 V for 30 min.

Step 11.

Block with 5% skim milk RT for 1 hour.

Step 12.

1:1000 Ab1 DnaA Ab.

Step 13.

4 °C over night rotation.

Step 14.

Wash the membrane with 1* PBST for 3 times (15 min , 10 min, 5 min).

Step 15.

Dillute AB2 1:10000 in 50 ml PBST.

Step 16.

Incubate with Ab2 for 1 h.

Step 17.

Wash the membrane with 1* PBST for 3 times (15 min , 10 min, 5 min).

Step 18.

Mix solution 1 and 2 (1:1).

Step 19.

Incubate the membrane with Mix solution for 5 s.

Step 20.

Develop the membrane with ChemiDoc MP Imaging system (BIO-RAD).
