

Aug 08,
2019

Primary Human Fibroblast Cell Culture

Aparna Kumar¹, Faria Zafar¹, Birgitt Schuele¹

¹Schuele Laboratory, Department of Pathology, Stanford University, Stanford, CA

1 Works for me dx.doi.org/10.17504/protocols.io.427gyhn

Aparna Kumar ⚡

ABSTRACT

Growing and maintaining primary human skin fibroblast cell cultures.

Additional links for protocols can be found:

- https://www.coriell.org/0/PDF/Fibroblast_Culture_FAQ.pdf
- <https://www.ncbi.nlm.nih.gov/pubmed/23852182>

GUIDELINES

Follow biosafety level 2 guidelines handling human tissue.

MATERIALS

NAME	CATALOG #	VENDOR
1.5ML Microcentrifuge Tube	BT620-NS	Bio Basic Inc.
EZ-LINE Serological Pipettes, 5ml stripette, Sterile, indiv wrapped, 250 per case	CR4113.SIZE.1CS	Bio Basic Inc.
DMSO	D8418	Sigma
Falcon® Serological Pipettes, 10 mL 200 Pipettes	38004	Stemcell Technologies
0.1% Gelatin in Water 500 mL	7903	Stemcell Technologies
Corning® Cryogenic Vials with Orange Caps, 2 mL 500 Vials	38053	Stemcell Technologies
FBS		Invitrogen - Thermo Fisher
Trypsin		Thermo Fisher Scientific
PBS		Invitrogen - Thermo Fisher
Trypan Blue Solution 0.4% (w/v) in PBS pH 7.5 ± 0.5	25-900-CI	Corning
AmnioMAX™ C-100 Basal Medium	17001074	Thermo Fisher Scientific
Remel™ 70% Ethanol	R2470110	Thermo Fisher Scientific
Nunc™ EasYFlask™ Cell Culture Flasks T25 Solid	156340	Thermo Fisher Scientific
Falcon™ 15mL Conical Centrifuge Tubes	14-959-53A	Fisher Scientific
Thermo Scientific™ Molecular BioProducts™ SoftFit-L Reload™ System: Sterile Tips	21-402-556	Fisher Scientific
Fisherbrand™ Aspirating Pipets	14-955-135	Fisher Scientific
Corning LX CoolCell Freezing System for Cryogenic Vials green 12 exposed vials; 1/ea	UX-04392-02	Cole Parmer
Fungizone (Amphotericin B)	15290018	Thermo Fisher Scientific
CHANG Medium C Frozen Supplement-14mL	C106	Fujifilm Wako Pure Chemical

Reagents

Medium, AmnioMax C100, kit
 0.1% Gelatin in Water, autoclaved
 Trypan Blue
 Ethanol 70%
 PBS
 FBS
 20% DMSO
 .25% Trypsin
 Fungizone (Amphotericin B)
 CHANG Medium C Frozen Supplement-14mL

Equipment

Labconco Biosafety Cabinet
 Benchtop Centrifuge
 Nanotek EVE Automated Cell Counter
 37 °C Water bath
 P1000 Rainin Pipet-Lite™
 P20 Rainin Pipet-Lite™
 -80 °C freezer
 Liquid Nitrogen tank

Materials

EVE Cell Counting Slide
 Falcon™ 15mL Conical Centrifuge Tubes
 50 mL Fisherbrand™ Sterile Polystyrene Disposable Serological Pipet
 25 mL Fisherbrand™ Sterile Polystyrene Disposable Serological Pipet
 10 mL Fisherbrand™ Sterile Polystyrene Disposable Serological Pipet
 5 mL Fisherbrand™ Sterile Polystyrene Disposable Serological Pipet
 Fisherbrand™ T25 flask
 Fisherbrand™ Premium Microcentrifuge Tubes: 1.5mL
 Thermo Scientific™ Molecular BioProducts™ SoftFit-L Reload™ System: Sterile Tips
 2 mL Fisherbrand™ Sterile Polystyrene Disposable Aspiration Pipet
 2mL Cryogenic Vials with Orange caps
 CoolCell Freezing System for Cryogenic Vials

SAFETY WARNINGS

Use appropriate safety measures when working with nitrogen tank and **-80 °C** freezer such as cryogenic gloves, full face visor and lab coats. Abide by biosafety level 2 standard safety procedures as human materials are used.

BEFORE STARTING







Turn on biosafety cabinet 10-15 min before start and clean all working surfaces with 70% ethanol. Make sure all materials in the cabinet are sterile and utilize standard aseptic technique.

Note: For plating density of human fibroblasts use $1.0-1.4 \times 10^4$ cells/cm².





Concepts and terms in Culturing Aging Fibroblasts

These include cumulative population doublings (CPD), senescence and senesced. CPD refers to the number of times that the cell number has been doubled. It has been used to measure the total number of cell divisions and it can be affected by several biological factors including the maximum lifespan of the species, age of the donor, the site of the biopsy, and the culture conditions. If cultures fail to reach confluency in 1 week, the culture is termed as "senescence". Cultures are considered to be "senescent" (at the end of their replicative lifespan) when they are unable to complete one population doubling during a 4-week period that includes 3 consecutive weeks of re-feeding with fresh medium containing 10% FBS. Reference: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3873382/>



Preparing the Media

- 1 Thaw  **15 ml** Chang C supplement and wipe bottle with 70% alcohol before placing in a biosafety cabinet.
- 2 Take the AmnioMax Basal Media and the Fungizone and wipe both bottles with 70% alcohol before placing in a biosafety cabinet.
- 3 Using a 50 ml pipette, aliquot  **90 ml** and using a 25 ml pipette, add the  **15 ml** Chang C supplement to the media.
- 4 Add  **100 µl** of Fungizone to the media.
- 5 Wrap the media container in foil to protect in the light and place in  **4 °C** until ready for use. Complete media can be stored at  **4 °C** for one month.





Preparing for Culturing

- 6 Store cryovials in liquid nitrogen storage tank upon arrival. Place vials in dry ice when ready to work with cells (do not thaw immediately).
 - 7 Record the label of the vial (ID, date frozen, any other documentation) and document in inventory.
 - 8 Prepare culture dish to thaw fibroblasts by adding  **2 ml** of 0.1% gelatin (autoclaved) to a culture dish and incubate at room temperature for  **00:20:00**.
- Note: For plating density of human fibroblasts use $1.0\text{--}1.4 \times 10^4$ cells/cm², e.g. for T25 flask plate 2.5×10^5 – 3.5×10^5 . If cells are plated too sparse the grow slows down and cells go into senescence.
- 9 Warm up AmnioMax media in water bath, wipe bottle with 70% alcohol before placing in a biosafety cabinet. Aliquot  **9 ml** of the AmnioMax media into a 15 ml conical tube.
 - 10 Aspirate gelatin and add  **3 ml** of AmnioMax media to the culture dish.



Culturing Cells

- 11 Thaw the cells by holding the lower half of the vial in  **37 °C** water bath and carefully swirling the vial around until only a small amount of ice is left in the vial. Do not let cells thaw completely. Wipe vial with 70% alcohol, before placing in biosafety cabinet.
- 12 Carefully remove the cap of vial and gently mix by triturating using a P1000 micropipette to resuspend the cells without causing foam.
- 13 Using a p20 micropipette, remove  **10 µl** from the vial for cell counting and add to a 1.5 ml microcentrifuge tube or onto a piece of

parafilm.



- 14 Using a p1000 micropipette, take  **400 µl** of media from the 15 ml conical tube and add to the vial with cells drop by drop. Note: depending on the size of the cryovial adjust the volume. Standard cryovials hold 1.5-2.0 ml.
- 15 Using a p1000 micropipette, aspirate solution from cryovial and slowly add to prepared conical tube with  **9 ml** Amniomax (Step 4) in a drop by drop fashion.
- 16 Centrifuge the 15 ml conical tube containing the cells and the medium at  **1000 rpm** for  **00:05:00** . Note: always balance centrifuge.

Cell counting

- 17 Take the 1.5 ml microcentrifuge tube or drop on parafilm from step 6 and add  **10 µl** of Trypan Blue and gently mix by pipetting up and down.
- 18 Take  **10 µl** of the solution and fill cell counting slide.
- 19 Insert the cell counting slide into the cell counter and count the cells. Record the number of cells and viability and then dispose of the cell counting slide.

Based on cell counts and viability adjust the culture vessel size to plate $1.0-1.4 \times 10^4$ of viable cells per cm^2 .

Culturing Cells

- 20 Take the 15 ml conical tube from the centrifuge and aspirate Amniomax media until there is only a little left above the cell pellet.
- 21 Add  **1 ml** of Amniomax media and gently triturate to resuspend the cell pellet.
- 22 Add 1ml cell solution from step 16 to the culture dish and gentle rock the culture dish back and forth to evenly distribute the cell over entire surface. Place flask in  **37 °C**, 5% CO₂ humidified incubator.

Maintaining Cell Culture

- 23 20 - 30% confluency should be reached 24 hours after the cells have been plated (the following day) 24hrs after plating. Figure Legend:
Example of a human fibroblast cell culture




Figure Legend: Example of a human fibroblast cell culture 24hrs after plating 1.2×10^4 cells per cm^2 .

- 24 Aspirate the media the cells are in and then add new media (the media should be warmed up and the bottles should be wiped down with 70% alcohol before use).
- 25 The whole media should be changed every 2 days until 90-100% confluency is reached, at which point the cells should be passaged.






Preparation for Passaging Cells

- 26 Add **4 ml** of 0.1% gelatin (autoclaved) to culture dishes (x2 T25) and incubate at room temperature for **00:20:00**. Note: For plating density of human fibroblasts use $1.0\text{--}1.4 \times 10^4$ cells/ cm^2 , e.g. for T25 flask plate 2.5×10^5 – 3.5×10^5 . If cells are plated too sparse the grow slows down and cells go into senescence.
- 27 Warm up AmnioMax media in water bath, wipe bottle with 70% alcohol before placing in a biosafety cabinet.
- 28 Aspirate gelatin and add **3 ml** of AmnioMax media to the culture dish.







Passaging the Cells

- 29 Remove Fibroblasts from incubator, place culture vessel in biosafety cabinet and aspirate media.
- 30 Add **2 ml** of PBS using 5ml serological pipette, swirl culture vessel, aspirate and discard PBS.
- 31 Add **2 ml** of 0.25% Trypsin and incubate for **00:05:00** at **37 °C** and 5% CO_2 . Firmly adherent cells can be detached quickly at 37°C . Observe the cells under microscope. The detached cells appear rounded under microscope. Tap the culture vessel gently to detach all cells.

Note: Prevent cell exposure to trypsin solution for longer periods (≥ 10 min)

- 32 Add  **4 ml** of media to neutralize the trypsin. Note: add 2 volumes of pre-warmed complete growth media (containing serum) to inactivate trypsin.
- 33 Using a pipette, completely mix cells with no clumps, transfer the cell suspension to a 15 ml conical centrifuge tube.
- 34 For cell counting, remove  **10 μ l** of solution from the tube using a P20 micropipette and add to a 1.5 ml microcentrifuge tube, see steps under "Cell Counting".
- 35 Centrifuge the cells in the 15 ml conical tube for  **00:05:00** at  **1000 rpm**.
- 36 Aspirate the supernatant until there is only a little left above the cell pellet.
- 37 Resuspend the cell pellet in media by gently triturating.
- 38 Remove 2.5×10^5 – 3.5×10^5 cells using a pipette (the amount of solution to pipette out is based on the number of cells which was determined by counting the cells) and add to each of the gelatin coated T25 flasks. Note: For plating density of human fibroblasts use 1.0 – 1.4×10^4 cells/cm², e.g. for T25 flask plate 2.5×10^5 – 3.5×10^5 . If cells are plated too sparse the grow slows down and cells go into senescence.
- 39 Place cells in incubator at  **37 °C** and 5% CO₂ and follow protocol for "Maintaining Cell Culture."

Cryofreezing Cells

- 40 If there are cells left in the 15 ml conical centrifuge tube that need to be cryofrozen, place the tube in the centrifuge and centrifuge for  **00:05:00** at  **1000 rpm**.
- 41 Aspirate any supernatant until there is only a little left above the cell pellet and resuspend the pellet by tapping the side of the tube.
- 42 Add  **500 μ l** of pre-warmed complete growth media (containing serum) to the conical tube with the cells using a P1000 micropipette and gently triturate.
- 43 Add  **500 μ l** of 20% DMSO (diluted in complete growth media) to the conical tube with the cells and gently triturate with a P1000 micropipette.
- 44 Using the pipette, transfer the cells to a Cryogenic Vial and label the vial. Note: If the time it will take between filling the vials and placing them in  **-80 °C** will be quite long, keep the vials in ice as they are being filled and transfer them to the freezing container once you are at the  **-80 °C** freezer.

- 45 Place the Cryogenic Vial in a Cryogenic Vial freezing system container and place in -80°C . Note: Make sure all of the holes in the freezing container are filled so that the temperature spreads evenly.
- 46 Transfer the vials from the -80°C to a Liquid Nitrogen tank no earlier than 24:00:00 after putting them in -80°C .

Results

47




Figure Legend: Fibroblasts stained with Fibronectin EP5 (555) by Santa Cruz Biotechnology (sc-8422) and Hoechst.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited