



Feb 27, 2019

Working

Thawing Rosettes and NPC

In 1 collection

Celeste Karch¹, Rita Martinez¹, Jacob Marsh¹¹Washington University in St Louis

dx.doi.org/10.17504/protocols.io.x9afr2e

Neurodegeneration Method Development Community

Tech. support email: ndcn-help@chanzuckerberg.com

Celeste Karch

Washington University in St Louis

IPSC CORTICAL
DIFFERENTIATION
022017.pdf

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

This protocol is part of the [IPSC CORTICAL DIFFERENTIATION](#) collection.

This method should be performed using sterile technique.


MATERIALS TEXT

Please refer to the attached full manuscript for required materials.

SAFETY WARNINGS

Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.

- 1 Aliquot 9 ml of pre-warmed DMEM/F12 into a 15 ml conical tube.
- 2 Remove cells from liquid nitrogen and thaw in 37 °C water bath for approximately 00:00:30 or until a small ice clump remains.
- 3 Add freshly thawed cells into the 9 ml of pre-warmed DMEM/F12 and mix gently by gently shaking tube 3 times. Avoid breaking up clumps of cells.
- 4 Centrifuge cells at 750 rpm for 00:03:00 . Aspirate supernatant.

5 Resuspend in  2 ml of NIM and plate into one well of a 6-well plate.



If cell viability count is necessary, after resuspension remove 10 μ l into a separate microfuge tube for cell viability count .



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited