



Jun 25,  
2019

## in situ hybridization performed on the peri-rhopalial tissue of a scyphozoan jellyfish [↗](#)

PLOS One

Christelle Bouchard<sup>1</sup>

<sup>1</sup>University of South Florida

Working

[dx.doi.org/10.17504/protocols.io.3qngmve](https://doi.org/10.17504/protocols.io.3qngmve)



Christelle Bouchard

### EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0218806>

### THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Bouchard C, Boudko DY, Jiang RHY (2019) A SLC6 transporter cloned from the lion's mane jellyfish (Cnidaria, Scyphozoa) is expressed in neurons. PLoS ONE 14(6): e0218806. doi: [10.1371/journal.pone.0218806](https://doi.org/10.1371/journal.pone.0218806)

- 1 Fix tissue 2 h, RT, in freshly prepared 4% PF in 0.1 M PB (pH 9.5) + 0.42M NaCl + 2 mM MgSO<sub>4</sub>
  - 1.1 Phosphate Buffer pH 9-9.5: made from phosphate dibasic (No pH adjustment required )
  - 1.2 PBS: 100 mM pH 7.5, 150 mM NaCl
- 2 Wash 3X 5' in PBST, pH 7.5
  - 2.1 PBST: PBS +0.1% Tween
- 3 Dehydration in PBS/methanol 3:1 → 1:1 → 1:3 10 min each then 100% methanol 5 min to ON in -20C.
- 4 Rehydration in PBS/methanol 1:3 → 1:1 → 3:1 10 min each then PBST, 10 min.
- 5 Heat treatment in 1 ml PBST, 5 min. Bain-marie at 95C with stopper. First 2 min, tissue should be resuspended vigorously every 20 sec, every 60 sec during the rest of incubation in order to prevent glueing of specimens. Snap-cool tubes on ice.
- 6 Wash 2X in PBST

- 7 Incubation 5 min in TEAA
  - 7.1 TEA: 0.1 M triethanolamine. Dissolve powder in water. Adjust pH at 7.8 with 5 N NaOH.  
TEAA: 5 ml TEA + 12, 5 ul glacial acetic acid. Acetic acid must be added to TEA just before the incubation with tissues.
- 8 Wash 2X 10 min in PBST
- 9 Post-fixation in 4% PF in PBST pH 7.5 for 20 min.
- 10 Wash 3X in PBST, 10 min each
- 11 Saturation of unspecific binding with 75 ul PBST/HS-DNA (0,5 mg/ml) for 10 min
  - 11.1 PBST/fish sperm: 0, 5 mg/ml fish sperm DNA in PBST. This DNA was made from herring and cod sperm and has a molecular weight distribution of 150-3000 bases.
- 12 Addition of the same amount of HybMix, incubate 5 min
  - 12.1 HybMix:  
50% formamide  
5 x SSC  
1mg/ml HS-DNA  
0, 1% Tween 20
- 13 Prehybridization in 150 ul HybMix for 10 min RT
- 14 Prehybridization in 150 ul HybMix for 2h at 45C
- 15 Hybridization in 1 ml HybMix + DIG-probes sense or antisense for 1 day and 2 nights at 40C. Before use, denature probe in HybMix at 70C for 10 min and snap-cool on ice.
- 16 2 wash in 1 ml of Wash I 60 min at 45C
  - 16.1 Wash I:  
50% formamide  
2 x SSC  
0, 1 % Tween 20
- 17 4 wash in 1 ml of Wash I, 15 min at 45C

- 18 Wash in PBST 20 min.
- 19 Saturation in PBST/10% NGS, 2h in fridge with shaking
- 19.1 Blocking solution: 10 % normal goat serum in PBST
- 20 Incubation in anti-DIG-AP-conjugated antibody diluted 1:2000 in PBST containing 1 % NGS/ 4C over the weekend with shaking. (preadsorption was not necessary).
- 21 Wash 4X in PBST for 20 min.
- 22 Incubation 3X 10min in AP-Buffer (containing 10mM MgCl<sub>2</sub>)
- 22.1
- 22.2
- 23 Substrate reaction with NBT/BCIP solution
- 24 Stop reaction with PBST/50 mM EDTA. 3X 5min
- 24.1
- 24.2
- 24.3
- 24.4
- 24.5
- 25

26

27

27.1



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited