

***S. cerevisiae* Δpho84 complementation: metal toxicity**

David Milner

Abstract

Testing for metal toxicity as a result of addition of a putative phosphate transporter to a *S. cerevisiae* Δpho84 strain

Citation: David Milner *S. cerevisiae* Δpho84 complementation: metal toxicity. **protocols.io**

dx.doi.org/10.17504/protocols.io.re4d3gw

Published: 25 Jul 2018

Protocol

Step 1.

Grow 10 ml *S. cerevisiae* strains BY5676 p416 GPD, BY2959 p416 GPD and BY2959 pXXX [p416 GPD + gene of interest] in Sc-ura medium for 16 hours at 30°C, 200 rpm shaking. Grow each strain in triplicate (3 independently obtained transformants).

Step 2.

Centrifuge cells at 1,810 x g for 3 minutes, wash twice and re-suspend in water.

Step 3.

Dilute cells to an OD₆₀₀ of 1.0 and prepare a dilution series (down to 10⁻⁴).

Step 4.

Spot 10 µl of each dilution onto Sc-ura plates containing each metal to be assayed. Use a Sc-ura plate with no added metal as a positive control to ensure that all strains grow equally well

Metals to try (try concentrations either side of these suggested starting points):

Metal	Suggested starting concentration
Potassium chromate	300 µM
Cobalt chloride	1 mM
Zinc chloride	10 mM
Sodium arsenate	3 mM
Manganese chloride	10 mM
Copper chloride	1 mM

Step 5.

Incubate plates for 72 hours at 30°C before imaging on a G:BOX imaging system (Syngene) to identify compounds on which toxicity was restored by the addition of pXXX.

Use these results to optimise compound concentrations for a second run.
