



Apr 16, 2019

Working

Immunoprecipitation assay [↗](#)

PLOS One

Ning Zhang , Peng Gao¹¹China Agricultural University

dx.doi.org/10.17504/protocols.io.zgbf3sn

China Agricultural University

Tech. support email: annachuang@cau.edu.cn

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0215408>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Zhang N, Gao P, Yin B, Li J, Wu T, Kuang Y, Wu W, Li J (2019) Cathepsin L promotes secretory IgA response by participating in antigen presentation pathways during *Mycoplasma Hyopneumoniae* infection. PLoS ONE 14(4): e0215408. doi: [10.1371/journal.pone.0215408](https://doi.org/10.1371/journal.pone.0215408)

GUIDELINES

IP-1.jpg
 IP-2.jpg

MATERIALS

NAME ▼	CATALOG # ▼	VENDOR ▼
PBS		Invitrogen - Thermo Fisher
Protein A/G Agarose		Santa Cruz Biotechnology
NP-40		Sigma Aldrich
MHC class II antibody	K274.3G8	Abcam

- 1 To harvest cells in nondenaturing conditions, remove media and rinse cells with ice cold 1×PBS
- 2 Remove PBS and add 0.5 ml ice cold cell lysis buffer to each plate(10cm) and incubate on ice for 5 min
- 3 native protein: Scrape cells off the plate and transfer to microcentrifuge tubes. Keep on ice. Sonicate on ice 3 times for 5 sec each. [denatured condition: scrape cells off the plate and transfer to 2.5 ml tube, boil for 10 min, then add 4 volumes cell lysis buffer, this is the cell lysate.]
- 4 Native protein : Microcentrifuge for 10 min at 4 °C,14000g and transfer the supernatant to a new tube. The supernatant is the cell lysates, store at -80°C.
- 5 add primary antibody(at the appropriate dilution as recommended in the product datasheet) to 200ul cell lysate at 1mg/ml. incubate with gentle rocking overnight at 4°C.

- 6 Add protein A/G agarose beads (10-30ul of 50% bead slurry). Incubate with gentle rocking for 1-3h at 4 °C.
- 7 Microcentrifuge for 30 sec at 4°C. Wash pellet five times with 500ul cell lysis buffer. Keep on ice between washes.
- 8 Resuspend the pellet with 20ul SDS sample buffer. Vortex, then microcentrifuge for 30sec.
- 9 Heat the sample to 95-100°C for 5 min and microcentrifuge for 1min at 14000g.Or, do not heat.
- 10 Load the sample on a 4-20% gel for SDS-PAGE.
- 11 Analyze sample by western blotting.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited