

Protocol for washing and modifying QCM-D sensors Version 2

Victoria Linderberg

Abstract

This protocol describes how to wash QCM-D sensors before usage and how to modify their surface using a composition of symmetric (poly)ethylene glycol (PEG) thiols consisting of 99% dS-PEG and 1% dS-PEG-biotin.

Citation: Victoria Linderberg Protocol for washing and modifying QCM-D sensors. **protocols.io**

dx.doi.org/10.17504/protocols.io.ex6bfre

Published: 20 May 2016

Before start

Protocol

Step 1.

Prepare solutions of 2% SDS, 0.1% SDS and 0.1 M HCl.

Step 2.

Clean the gold-coated QCM-D sensors by hand with 2% SDS to mechanically remove contaminants. Be sure that the glove is wet all the time so it doesn't scratch the sensor!

Step 3.

Place the sensors in a stand and put them in a 0.1% SDS solution. Heat for 10 min at 80°C, then rinse the sensors with water.

Step 4.

Place the sensors in a 0.1 M HCl solution. Heat for 10 min at 80°C, then rinse the sensors with water several times to remove all HCl solution. Then blow the sensors dry under a flow of nitrogen.

Step 5.

Make an ethanolic solution with total concentration 10 mM of mixed PEG thiols (PEG and biotin-PEG thiols in a molar ratio of 99:1).

($M_{W, dS-PEG} = 386.5$ Da and $M_{W, dS-PEG-biotin} = 788.0$ Da)

Step 6.

Incubate the QCM-D sensors in the ethanolic PEG thiol solution for at least 16 hours at 4°C in darkness.

Step 7.

Rinse the sensors a couple of times with ethanol, then a couple of times with milli-Q water and lastly blow them dry with nitrogen. They may now be used or stored.

Warnings

If storing the modified sensors, be sure to keep them cold and in darkness!