

Lminex Milliplex Soluble Cytokine Receptor 13-plex

Troy Kemp, Ligia Pint

Abstract

Lminex Milliplex Soluble Cytokine Receptor 13-plex manufacturer's protocol

Citation: Troy Kemp, Ligia Pint Lminex Milliplex Soluble Cytokine Receptor 13-plex. **protocols.io**

dx.doi.org/10.17504/protocols.io.hvib64e

Published: 01 Aug 2017

Protocol

Step 1.

PREPARATION OF SAMPLES/REAGENTS FOR IMMUNOASSAY

Step 2.

Preparation of Serum/Plasma Thaw Time: Thaw the samples completely on ice, mix well by vortexing and centrifuge (10,000 rpm, 10 minutes, 4°C) prior to use in the assay to remove particulates. Serum/Plasma samples from normal subjects should be diluted 1:5 using the Serum Matrix provided in the kit as the sample diluent (20 µL sample mixed with 80 µL Serum Matrix). If samples require dilution beyond 1:5, continue to use the Serum Matrix as the sample diluent.

Step 3.

Preparation of Antibody-Immobilized Beads For individual vials of beads, sonicate each antibody-bead vial for 30 seconds then vortex for 1 minute. Add 60 µL from each antibody-bead vial to the Mixing Bottle and bring final volume to 3.0 mL with Bead Diluent. Vortex the mixed beads well. Unused portion may be stored at 2-8°C for up to one month. Example: When using 13 antibody-immobilized beads, add 60 µL from each of the 13 bead sets to the Mixing Bottle. Then add 2.22 mL Bead Diluent.

Step 4.

Preparation of Quality Controls Reconstitution Time: Before use, reconstitute Quality Control 1 and Quality Control 2 with 250 µL deionized water. Invert the vial several times to mix and vortex. Allow the vial to sit for 5-10 minutes and then transfer the controls to appropriately labeled polypropylene microfuge tubes. Unused portion may be stored at -20°C for up to one month.

Step 5.

Preparation of Wash Buffer Bring the 10X Wash Buffer to room temperature and mix to bring all salts into solution. Dilute 30 mL of 10X Wash Buffer with 270 mL deionized water. Store unused portion at 2-8°C for up to one month.

Step 6.

Preparation of Serum Matrix Reconstitution Time: Add 5.0 mL Assay Buffer to the bottle containing lyophilized Serum Matrix. Mix well. Allow at least 10 minutes for complete reconstitution. Leftover reconstituted Serum Matrix should be stored at -20°C for up to one month.

Step 7.

Preparation of Human Soluble Cytokine Receptor Standard Reconstitution Time: 1.) Prior to use, reconstitute the Human Soluble Cytokine Receptor Panel Standard with 250 µL deionized water to

Add 25 µL of the Serum Matrix solution to the background, appropriate standards, and control wells.

Step 19.

Add 25 µL of Sample into the appropriate wells.

Step 20.

Vortex Mixing Bottle and add 25 µL of the mixed Beads to each well. (Note: During addition of Beads, shake bead bottle intermittently to avoid settling.)

Step 21.

Seal the plate with a plate sealer, cover it with the lid. Wrap a rubber band around the plate holder, plate and lid and incubate with agitation on a plate shaker 2 hours at room temperature (20-25°C).

Step 22.

Gently remove fluid by vacuum.

Step 23.

Wash plate 2 times with 200 µL/well of Wash Buffer, removing Wash Buffer by vacuum filtration between each wash. Blot excess Wash Buffer from the bottom of the plate with an absorbent pad or paper towels.

Step 24.

Add 25 µL of Detection Antibodies into each well. (Note: Allow the Detection Antibodies to warm to room temperature prior to addition.)

Step 25.

Seal, cover with lid, and incubate with agitation on a plate shaker for 1 hour. DO NOT VACUUM AFTER INCUBATION.20) Add 25 µL Streptavidin-Phycoerythrin to each well containing the 25 µL of Detection Antibodies.

Step 26.

Seal, cover with lid and incubate with agitation on a plate shaker for 30 minutes at room temperature (20-25°C).

Step 27.

Gently remove all contents by vacuum.

Step 28.

Wash plate 2 times with 200 µL/well Wash Buffer, removing Wash Buffer by vacuum filtration between each wash. Wipe any excess buffer on the bottom of the plate with a tissue.

Step 29.

Add 150 µL of Sheath Fluid to all wells. Resuspend the beads on a plate shaker for 5 minutes.

Step 30.

Run plate on Luminex 100™ IS.

Step 31.

Save and analyze the data using Bio-Plex Manager software.EQUIPMENT SETTINGSEvents: 50, per bead region Sample Size: 100 µLGate Settings 4335 to 10,000Time Out 60 secondsQUALITY CONTROLThe ranges for each analyte in Quality Control 1 and 2 are provided on the card insert or can be located at the MILLIPORE website www.millipore.com/techlibrary/index.do using the catalog number as the keyword.