

# Q5® Site-Directed Mutagenesis (E0552) Version 2

# **New England Biolabs**

# **Abstract**

This is the protocol for the Q5® Site-Directed Mutagenesis Kit without competent cells

Citation: New England Biolabs Q5® Site-Directed Mutagenesis (E0552). protocols.io

dx.doi.org/10.17504/protocols.io.is3cegn

Published: 02 Mar 2018

# **Materials**

Q5 Site-Directed Mutagenesis Kit (Without Competent Cells) - 10 rxns <u>E0552S</u> by <u>New England</u> Biolabs

# **Protocol**

# **Exponential Amplification (PCR)**

#### Step 1.

Assemble the following reagents in a thin-walled PCR tube.

12 E	
12.5 μι	1X
1.25 μl	0.5 μΜ
1.25 μl	0.5 μΜ
1 μΙ	1-25 ng
9.0 μΙ	
	1.25 μl 1 μl



# . E0552 Q5 PCR Mixture

**CONTACT: New England Biolabs** 

Step 1.1.

Q5 Hot Start High-Fidelity 2X Master Mix 12.5 μl

# Step 1.2.

10 μM Forward Primer 1.25μl

# Step 1.3.

10 μM Reverse Primer 1.25μl

# Step 1.4.

Template DNA (1-25 ng/μl) 1μl

# Step 1.5.

Nuclease-free water 9µl

# Exponential Amplification (PCR)

# Step 2.

Mix reagents completely.

# Exponential Amplification (PCR)

## Step 3.

Transfer to a thermocycler and perform the following cycling conditions:

Thermocycling Conditions for a Routine PCR:

STEP	TEMP	TIME	
Initial Denaturation	98°C	30 seconds	
	98°C	10 seconds	
25 Cycles	50-72°C*	10-30 seconds	
•	72°C	20-30 seconds/kb	
Final Extension	72°C	2 minutes	
Hold	4-10°C		

#### NOTES

# New England Biolabs 09 Oct 2014

# Kinase, Ligase & DpnI (KLD) Treatment

# Step 4.

Assemble the following reagents:

	VOLUME FINAL CONC.		
PCR Product	1 μΙ		
2X KLD Reaction Buffer	r 5 μl	1X	
10X KLD Enzyme Mix	1 μΙ	1X	
Nuclease-free Water	3 μΙ		



<sup>\*</sup> For a Q5-optimized annealing temperature of mutagenic primers, please use  $\underline{\text{NEBaseChanger}^{\text{TM}}}$ , the online NEB primer design software. For pre-designed, back-to-back primer sets, a Ta = Tm + 3 rule can be applied, but optimization may be necessary.

#### . E0552 KLD Mixure

**CONTACT: New England Biolabs** 

Step 4.1.

PCR Product 1µl

Step 4.2.

2X KLD Reaction Buffer 5µl

Step 4.3.

10X KLD Enzyme Mix 1μl

Step 4.4.

Nuclease-free Water 3µl

# Kinase, Ligase & DpnI (KLD) Treatment

#### Step 5.

Mix well by pipetting up and down.

# Kinase, Ligase & DpnI (KLD) Treatment

#### Step 6.

Incubate at room temperature for 5 minutes.

© DURATION

00:05:00

# Transformation

#### Step 7.

Thaw a  $50~\mu l$  aliquot of chemically competent E. coli cells on ice.

#### NOTES

# New England Biolabs 26 Jan 2015

NEB 5-alpha Competent E. coli (High Efficiency), NEB #C2987, are recommended

# **Transformation**

# Step 8.

Add 5 µl of the KLD mix from the "KLD Section" above to the tube of thawed cells.

# **Transformation**

## Step 9.

Carefully flick the tube 4-5 times to mix. Do not vortex.

## **Transformation**

#### Step 10.

Place the mixture on ice for 30 minutes.

© DURATION

00:30:00

# **Transformation**

#### **Step 11.**

Heat shock at 42°C for 30 seconds.

**O DURATION** 

00:00:30

# **Transformation**

# Step 12.

Place on ice for 5 minutes.

**O DURATION** 

00:05:00

# Transformation

#### **Step 13.**

Pipette 950  $\mu$ l of room temperature SOC into the mixture.

AMOUNT

950 µl Additional info:



**REAGENTS** 

SOC Outgrowth Medium - 100 ml <u>B9020S</u> by New England Biolabs

## **Transformation**

# **Step 14.**

Incubate at 37°C for 60 minutes with shaking (250 rpm).

**O** DURATION

01:00:00

# **Transformation**

#### **Step 15.**

Mix the cells thoroughly by flicking the tube and inverting.

#### **Transformation**

#### **Step 16.**

Spread 50-100 µl onto a selection plate.

# **Transformation**

# **Step 17.**

Incubate overnight at 37°C

O DURATION

15:00:00

#### NOTES

## New England Biolabs 03 Oct 2014

It may be necessary (particularly for simple substitution and deletion experiments) to make a 10to 40-fold dilution of the transformation mix in SOC prior to plating, to avoid a lawn of colonies