



Oct 16, 2019

Removal of gDNA from totalRNA

[iGEM Dusseldorf¹](#)¹Heinrich-Heine Universität Düsseldorf**1** *Works for me* [dx.doi.org/10.17504/protocols.io.8bfhsjn](https://doi.org/10.17504/protocols.io.8bfhsjn) [iGEM Dusseldorf](#) 

ABSTRACT

Removal of genomic DNA from RNA preparations

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

https://assets.thermofisher.com/TFS-Assets/LSG/manuals/MAN0012000_DNase_I_RNasefree_1UuL_UG.pdf

GUIDELINES

Always work with gloves and safety gear and work on ice.

MATERIALS TEXT


- RNase-free H₂O
- DNaseI-buffer (10X)
- DNaseI
- 50 mM EDTA
- NaOAc (3M & pH 5,3) and 100% EtOH

Digestion of DNA

- 1 Calculate the volume required for 1 µg RNA

Component	amount
DNaseI buffer	1 µl
total RNA	1 µg
DNaseI	1 µl
H2O	to 10 µl

- 2 pipet components into a RNase-free 1.5 ml tube (following order: H₂O, buffer, total RNA, DNaseI)
- 3 incubate reaction at 37°C

 **00:30:00**

Extraction of DNaseI-digested RNA

- 4 add 1 µl 50 mM EDTA and incubate at 65°C

🕒 00:10:00

- 5 add 1/10 Volume NaOAc and 3 Volumes 100% EtOH

- 6 freeze at -20°C over night or at -80°C for at least 1h

- 7 centrifuge at 13000 rpm at 4°C

🕒 00:30:00

- 8 Discard supernatant

RNA pellet may be invisible, make sure to align tube lid so you remember where the pellet is located

- 9 add 500 µl 70% EtOH

- 10 centrifuge at 13000 rpm at 4°C

🕒 00:10:00

- 11 discard supernatant, remove as much ethanol as possible

- 12 dry tubes at RT under the hood (~5min), do not overdry!

- 13 resuspend pellet in 20 µl RNase-free water

- 14 measure concentration

Quality control:

- 15 perform PCR on your DNA-free RNA. Include a positive control using gDNA as a template and primer that bind to it.

There should be no visible bands except for positive control!



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited