

# Fixation of marine samples for flow cytometry sorting

Daniel Vaultot

## Abstract

Protocol to fix marine samples for flow cytometry sorting of phytoplankton.  
Fix at least 2 samples per depth sampled and at least 6 to 10 depths per vertical profile.

Back to the laboratory, pico and nano-plankton populations can be sorted and used for clone library construction or metabarcoding with Next Generation Sequencing (e.g. Illumina). This was successfully tested on samples from a cruise off Brazil (unpublished data)

## Reference

Ribeiro C., Lopes A., Marie D., Vaultot D. Unpublished data

**Citation:** Daniel Vaultot Fixation of marine samples for flow cytometry sorting. **protocols.io**  
dx.doi.org/10.17504/protocols.io.d2x8fm

**Published:** 20 Oct 2015

## Before start

### Necessary equipment

- Vortex mixer
- Cryotubes 2 mL
- Liquid nitrogen tank

## Materials

Pluronic 10% [P5556-100mL](#) by [Sigma Aldrich](#)

DMSO [472301](#) by [Sigma Aldrich](#)

## Protocol

### Step 1.

Prefilter seawater sample onto 200 µm mesh

 **AMOUNT**

2 ml Additional info:

### Step 2.

Add 1.5 mL of prefiltered seawater sample to a 2 mL cryotube

### Step 3.

Add DMSO

 **AMOUNT**

150 µl Additional info:



#### REAGENTS

DMSO [472301](#) by [Sigma Aldrich](#)

#### Step 4.

Add Pluronic (facultative)



#### AMOUNT

2 µl Additional info:



#### REAGENTS

Pluronic 10% [P5556-100mL](#) by [Sigma Aldrich](#)

#### Step 5.

Vortex

#### Step 6.

Wait 10 min



#### DURATION

00:10:00

#### Step 7.

Flash freeze in liquid nitrogen tank

## Warnings

Samples must be stored either in liquid nitrogen or at -80°C, not at -20°C because degradation will take place at the latter temperature