SYBR Gold Staining for Viral Enumeration (Case 2)

Li Deng

Abstract

Case 2: Fix samples. Use this protocol when long term storage of slide is required.

For when you can count your samples in a few days, see <u>Case 1</u>.

Citation: Li Deng SYBR Gold Staining for Viral Enumeration (Case 2). protocols.io

dx.doi.org/10.17504/protocols.io.c7fzjm

Published: 04 Jan 2016

Protocol

Step 1.

Prepare buffer glutaraldehyde:



Buffered glutaraldehyde

CONTACT: Bonnie Poulos

Step 1.1.

Prepare solution A and solution B as found in guidelines

Step 1.2.

Prepare the phosphate buffere by adding 15.25 ml solution A and 9.75 ml solution B in 25 ml 0.02 μm filtered mQ

NOTES

Bonnie Poulos 24 Jun 2015

= 50 ml buffer

Step 1.3.

Add 40 ml of 25% glutaraldehyde to 10 ml phosphate buffer for a 20% buffered glutaraldehyde solution

NOTES

Bonnie Poulos 30 Nov 2015

This is most often used for transmission electron microscopy so purchase a TEM-grade glutaraldehyde (from a vendor such as Ted Pella) for preparing this buffered reagent

Step 1.4.

Adjust pH to 7 with NaOH

Step 1.5.

Store this 20% buffered glutaraldehyde at 4 °C in the dark

Step 2.

Make dilution of virus prep in 0.02 µm filtered seawater to a concentration of E+07 particles ml⁻¹

Step 3.

Fix viral sample with 20% buffered glut to a final concentration of 0.5%

Step 4

Incubate at 4°C for 30 min.

© DURATION

00:30:00

Step 5.

Prepare working solutions of SYBR Gold

PROTOCOL

. **SYBR Gold working solutions**

CONTACT: VERVE Team

Step 5.1.

Thaw the commercial stock of SYBR Gold in the dark at RT

Step 5.2.

Centrifuge at 3000 rpm for 5 minutes

O DURATION

00:05:00

NOTES

VERVE Team 24 Jun 2015

Because SYBR Gold is in DMSO.

Step 5.3.

Dilute SYBR Gold in 0.02 μm filtered TE buffer to 100x (10 μl in 990 μl TE buffer)

Step 5.4.

0.02 µm filter the diluted SYBR (100x)

NOTES

VERVE Team 24 Jun 2015

This working stock can be stored at -20°C and re-thawed one time.

Step 6.

Add 2 µl of SYBR working stock in 98 µl 0.02 µm filtered mQ in a plastic Petri dish, 4 drops in one dish

Step 7.

Cover the dish by aluminum foil

NOTES

VERVE Team 24 Jun 2015

SYBR is light sensitive.

Step 8.

Set up the filtration unit, connecting it to a vacuum

NOTES

VERVE Team 24 Jun 2015

Set up the vacuum no higher than 5 mm Hg.

Step 9.

Add a few drops of 0.02 µm filtered mQ on the filter base

Step 10.

Place a 0.2 nitrocellulose filter (the support filter) on top of the water

Step 11.

Switch on the vacuum, the support filter should be flat on the filter base

NOTES

VERVE Team 11 Aug 2015

This support filter is good for several samples as long as it remains flat and no air bubbles between filter and base.

Step 12.

Add a few drops of 0.02 µm filtered mQ on the support filter

Step 13.

Switch on the vacuum to pull the water through

Step 14.

Apply a 0.02 µm Anodisc filter over the support filter

NOTES

VERVE Team 11 Aug 2015

Make sure no air bubbles between filters.

Step 15.

Apply the filter tower and clamp while vacuum is on

Step 16.

Switch off the vacuum and add viral samples

Step 17.

Switch on the vacuum and wash filter set with 1 ml of 0.02 µm filtered mQ

■ AMOUNT

1 ml Additional info:

Step 18.

Remove the filter while the vacuum is still on

Step 19.

Rinse tower in 1L Q-water in between samples

Step 20.

Blot onto paper towel to dry

Step 21.

Dry filter membrane on Kimwipes at RT completely

Step 22.

Remove membrane and place viruses side up on staining solution in the Petri dish for 15 minutes (cover with aluminum foil)

© DURATION

00:15:00

Step 23.

Dry filter membrane again on Kimwipes in the **dark** at RT completely

NOTES

VERVE Team 26 Jun 2015

Better in a paper box.

Step 24.

Pipet 20 µl antifade solution on a microscope slide

Step 25.

Place the stained filter membrane on top of it

Step 26.

Pipet 30 µl antifade solution on a cover slide

Step 27.

Carefully place it on the filter to avoid bubbles

Step 28.

Place slide at -20°C to enhance fluorescence

Step 29.

Read slides using 100x oil immersion objective and inverted fluorescent microscope