



Jun 12,
2019

Rice medium

Roscoff Culture Collection¹

¹CNRS-Sorbonne Université, Station Biologique, Place G. Tessier 29680 Roscoff FRANCE

Working

[dx.doi.org/10.17504/protocols.io.s35egq6](https://doi.org/10.17504/protocols.io.s35egq6)

Roscoff Culture Collection



Daniel Vaultot

Station Biologique, Roscoff, France



ABSTRACT

Used for heterotrophic protists.

- 1 Filter 1L of seawater aged for at least two months on prefilter and filter 0.2 μ m
- 2 Heat seawater during 20 min at 100°C
- 3 Add 40 graism of rice
- 4 Autoclave 3 times the medium (let temperature decrease between each round of autoclaving)



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited