

# Luminex Milliplex Cytokine/Chemokine 9-plex MAG

Ligia Pinto, Troy Kemp

## Abstract

Luminex Milliplex Cytokine/Chemokine 9-plex MAG manufacturer's protocol

**Citation:** Ligia Pinto, Troy Kemp Luminex Milliplex Cytokine/Chemokine 9-plex MAG. [protocols.io](https://doi.org/10.17504/protocols.io.hvvhb636)

[dx.doi.org/10.17504/protocols.io.hvvhb636](https://doi.org/10.17504/protocols.io.hvvhb636)

**Published:** 01 Aug 2017

## Protocol

### Step 1.

PREPARATION OF SAMPLES/REAGENTS FOR IMMUNOASSAY

### Step 2.

Preparation of Serum/Plasma Thaw Time: Thaw the samples completely on ice, mix well by shaking on plate shaker for 1 min. at RT (20-25°C) and centrifuge (1,700 xg, 10 minutes, 4°C) prior to use in the assay to remove particulates.

### Step 3.

Preparation of Antibody-Immobilized Beads Sonicate each antibody-bead vial for 30 seconds; vortex for 1 minute. Add 150 µL from each antibody bead vial to the Mixing Bottle and bring final volume to 3.0 mL with Bead Diluent. Vortex the mixed beads well. Unused portions may be stored at 2-8°C for up to one month. Example: When using 9 antibody-immobilized beads, add 150 µL from each of the 9 bead sets to the Mixing Bottle. Then add 1.65 mL Bead Diluent.

### Step 4.

Preparation of Quality Controls      Reconstitution Time: Before use, reconstitute Quality Control 1 and Quality Control 2 with 250 µL deionized water. Invert the vial several times to mix and vortex. Allow the vial to sit for 5-10 minutes and then transfer the controls to appropriately labeled polypropylene microfuge tubes. Unused portion may be stored at £ -20°C for up to one month.

### Step 5.

Preparation of Wash Buffer Bring the 10X Wash Buffer to room temperature and mix to bring all salts into solution. Dilute 30 mL of 10X Wash Buffer with 270 mL deionized water. Store unused portion at 2-8°C for up to one month.

### Step 6.

Preparation of Serum Matrix      Reconstitution Time: Add 1.0 mL Assay Buffer to the bottle containing lyophilized Serum Matrix. Mix well. Allow at least 10 minutes for complete reconstitution. Leftover reconstituted Serum Matrix should be stored at £ -20°C for up to one month.

### Step 7.

Preparation of Human Metabolic Hormone Panel Standard Reconstitution Time: 1.) Prior to use, reconstitute the Human Metabolic Hormone Panel Standard with 250 µL deionized water to give STD7. Invert the vial several times to mix. Vortex the vial for 10 seconds. Allow the vial to sit for 5-10 minutes and then transfer the standard to an appropriately labeled polypropylene microfuge tube.

This standard will be termed STD7; the unused portion may be stored at  $-20^{\circ}\text{C}$  for up to one month.

#### **Step 8.**

) Preparation of Working Standards. Label six polypropylene microfuge tubes STD6, STD5, STD4, STD3, STD2, and STD1. Add 200  $\mu\text{L}$  of Assay Buffer to each of the six tubes. Prepare serial dilutions by adding 100  $\mu\text{L}$  of STD7 reconstituted standard to the STD6 tube, mix well and transfer 100  $\mu\text{L}$  of the STD6 standard to the STD5 tube, mix well and transfer 100  $\mu\text{L}$  of the STD5 standard to the STD4 tube, mix well and transfer 100  $\mu\text{L}$  of the STD4 standard to STD3 tube, mix well and transfer 100  $\mu\text{L}$  of the STD3 standard to the STD2 tube and mix well, transfer 100  $\mu\text{L}$  of the STD2 standard to the STD1 tube and mix well. The 0 pg/mL standard (Background) will be Assay Buffer. Standard Volume of Deionized Water to Add (mL) Volume of Standard to Add Original (STD7) 2500 Standard Concentration (pg/mL) Volume of Assay Buffer to Add (mL) Volume of Standard to Add STD6 200 100 mL of STD7 STD5 200 100 mL of STD6 STD4 200 100 mL of STD5 STD3 200 100 mL of STD4 STD2 200 100 mL of STD3 STD1 200 100 mL of STD2 After dilution, each tube has the following concentrations for each analyte: Standard Tube # GIP (pg/mL) Ghrelin, GLP-1, Glucagon, PP, PYY (pg/mL) Amylin (pg/mL) C-Peptide (pg/mL) Insulin, Leptin (pg/mL) 1

#### **Step 9.**

7

#### **Step 10.**

7

#### **Step 11.**

4

#### **Step 12.**

6

#### **Step 13.**

22

#### **Step 14.**

2

#### **Step 15.**

2

#### **Step 16.**

3205.8

#### **Step 17.**

53

#### **Step 18.**

7

#### **Step 19.**

5

#### **Step 20.**

9

#### **Step 21.**

31,2354

#### **Step 22.**

1370.4740.71,8523,7045

#### **Step 23.**

21,1112,2225,55611,1116

**Step 24.**

73,3336,66716,66733,33372,00010,00020,00050,000100,000IMMUNOASSAY PROCEDURE□ Allow all reagents to warm to room temperature (20-25°C) before use in the assay.□ Run the standards, controls, and samples in duplicate.

**Step 25.**

Prewet plate by pipetting 200 µL of Assay Buffer into each well of the MAG Plate. Seal and mix on a plate shaker for 10 minutes at room temperature (20-25°C).

**Step 26.**

Decant Assay Buffer and remove residual amount from all wells by inverting the plate and tapping it smartly onto absorbent towels several times.10) Add 25 µL of each Standard or Control into the appropriate wells. Add 25 µL Assay Buffer to the 0 pg/mL standard (Background).

**Step 27.**

Add 25 µL of Assay Buffer to the sample wells.

**Step 28.**

Add 25 µL of the Serum Matrix solution to the background, appropriate standards, and control wells.

**Step 29.**

Add 25 µL of Sample into the appropriate wells.

**Step 30.**

Vortex Mixing Bottle and add 25 µL of the mixed Beads to each well. (Note: During addition of Beads, shake bead bottle intermittently to avoid settling. Due to the composition of magnetic beads, you may notice a slight color in the bead solution. This does not affect the performance of the beads or the kit).

**Step 31.**

Seal the plate with a plate sealer, cover it with the lid. Wrap a rubber band around the plate holder, plate and lid and incubate with agitation on a plate shaker overnight (16-18 hours) at 4°C.

**Step 32.**

Gently remove fluid by aspiration.

**Step 33.**

Wash plate 3 times with 200 µL/well of Wash Buffer, removing Wash Buffer by aspiration between each wash.

**Step 34.**

Add 50 µL of Detection Antibodies into each well. (Note: Allow the Detection Antibodies to warm to room temperature prior to addition.)

**Step 35.**

Seal, cover with lid, and incubate with agitation on a plate shaker for 30 minutes at room temperature. DO NOT WASH AFTER INCUBATION.20) Add 50 µL Streptavidin-Phycoerythrin to each well containing the 50 µL of Detection Antibodies.

**Step 36.**

Seal, cover with lid and incubate with agitation on a plate shaker for 30 minutes at room temperature (20-25°C).

**Step 37.**

Gently remove all contents by aspiration.

**Step 38.**

Wash plate 3 times with 200 µL/well Wash Buffer, removing Wash Buffer by aspiration between each wash.

**Step 39.**

Add 100 µL of Sheath Fluid to all wells. Resuspend the beads on a plate shaker for 5 minutes.

**Step 40.**

Run plate on Luminex 100™ IS.

**Step 41.**

Save and analyze the data using Bio-Plex Manager software. EQUIPMENT SETTINGS  
Events: 50, per bead region  
Sample Size: 50 µL  
Gate Settings 5000 to 25,000  
Time Out 60 seconds  
QUALITY CONTROL  
The ranges for each analyte in Quality Control 1 and 2 are provided on the card insert or can be located at the MILLIPORE website [www.millipore.com/techlibrary/index.do](http://www.millipore.com/techlibrary/index.do) using the catalog number as the keyword.