

# Extraction method A

Sarah Siu Tze Mak, Shyam Gopalakrishnan, Christian Caroe, Chunyu Geng, Shanlin Liu, Mikkel-Holger S Sinding, Lukas F K Kuderna, Wenwei Zhang, Shujin Fu, Filipe G Vieira, Mietje Germonpré, Hervé Bocherens, Sergey Fedorov, Bent Petersen, Thomas Sicheritz-Ponten, Tomas Marques-Bonet, Guojie Zhang, Hui Jiang, M Thomas P Gilbert

## Abstract

This protocol provides an efficient DNA extraction and purification of historical museum hides, which potentially have been chemically tanned.

**Citation:** Sarah Siu Tze Mak, Shyam Gopalakrishnan, Christian Caroe, Chunyu Geng, Shanlin Liu, Mikkel-Holger S Sinding, Lukas F K Kuderna, Wenwei Zhang, Shujin Fu, Filipe G Vieira, Mietje Germonpré, Hervé Bocherens, Sergey Fedorov, Bent Petersen, Thomas Sicheritz-Ponten, Tomas Marques-Bonet, Guojie Zhang, Hui Jiang, M Thomas P Gilbert Extraction method A. **protocols.io**

dx.doi.org/10.17504/protocols.io.iabcaan

**Published:** 14 Jun 2017

## Before start

Clean

## Materials

- ✓ proteinase K containing buffer by Contributed by users
- ✓ organic solvents (phenol:chloroform) by Contributed by users
- MinElute columns by [Qiagen](#)
- PB buffer [19066](#) by [Qiagen](#)
- Buffer PE [19065](#) by [Qiagen](#)
- Buffer EB [19086](#) by [Qiagen](#)

## Protocol

### Extraction

#### Step 1.

Digestion in a proteinase K containing buffer following Gilbert et al. 2007.

## NOTES

**GigaScience Database** 02 Jun 2017

Gilbert MTP, Tomsho LP, Rendulic S, Packard M, Drautz DI, Sher A, et al. Whole-genome shotgun sequencing of mitochondria from ancient hair shafts. Science. American Association for the Advancement of Science; 2007;317:1927–30.

## Extraction

### Step 2.

Pre-digest all samples at 56 °C for 1 hour.

 **DURATION**

01:00:00

 **NOTES**

**GigaScience Database** 02 Jun 2017

Perform incubation with a rotation of the samples in an oven.

## Extraction

### Step 3.

Replace the buffer with fresh buffer.

## Extraction

### Step 4.

Perform a second 12-hour digest with the fresh buffer.

 **DURATION**

12:00:00

 **NOTES**

**GigaScience Database** 02 Jun 2017

Perform incubation with a rotation of the samples in an oven.

## Extraction

### Step 5.

Purify digests using organic solvents (phenol:chloroform) and Qiagen MinElute columns (Qiagen, Hilden, DE) as described in the following steps:

## Extraction

### Step 6.

Mix digests with 8x modified PB buffer (Allentoft et al. 2015).

 **NOTES**

**GigaScience Database** 02 Jun 2017

Allentoft ME, Sikora M, Sjögren K-G, Rasmussen S, Rasmussen M, Stenderup J, et al. Population genomics of Bronze Age Eurasia. Nature. Nature Publishing Group; Dabney J, Knapp M, Glocke I et al. (2013) Complete mitochondrial genome sequence of a Middle Pleistocene cave bear reconstructed from ultrashort DNA fragments. Proceedings of the National Academy of Sciences of the United States of America, 110, 15758–15763.2015;522:167–72.

#### Extraction

##### Step 7.

Bind the digests to the MinElute column using reservoirs as previously described (Dabney et al. 2013).

#### Extraction

##### Step 8.

Wash with 800 µl buffer PE (Qiagen).

#### Extraction

##### Step 9.

Elute in 17 µl buffer EB (Qiagen). (1/2)

#### Extraction

##### Step 10.

Incubate for 5 minutes at 37 °C prior to release. (1/2)

 DURATION

00:05:00

#### Extraction

##### Step 11.

Elute in 17 µl buffer EB (Qiagen). (2/2)

#### Extraction

##### Step 12.

Incubate for 5 minutes at 37 °C prior to release. (2/2)

 DURATION

00:05:00

#### Extraction

##### Step 13.

Prior to library construction, analyze small aliquots of each extract on an Agilent 2200 TapeStation HS chip (Agilent Technologies, Palo Alto, California, USA) for fragment size estimation and molar concentration.