

Extraction of genomic DNA from diatoms by the modified method described in (Jacobs et al., 1992)

Artyom M. Marchenkov, Darya P. Petrova, Alexey A. Morozov, Yulia R. Zakharova, Michael A. Grachev, Alexander A. Bondar

Abstract

A protocol for diatom DNA extraction used in the Limnological Institute SB RAS for freshwater diatoms.

Citation: Artyom M. Marchenkov, Darya P. Petrova, Alexey A. Morozov, Yulia R. Zakharova, Michael A. Grachev, Alexander A. Bondar Extraction of genomic DNA from diatoms by the modified method described in (Jacobs et al., 1992). **protocols.io**

dx.doi.org/10.17504/protocols.io.qh6dt9e

Published: 09 Jun 2018

Protocol

Cells collecting

Step 1.

The cells were sedimented on polycarbonate filters (5- μ m pores) (Whatman, USA), briefly rinsed with cold Diatom Medium (DM), harvested by centrifugation for 2 min at 16,100 g at 4 °C and then stored at -70 °C.

Cells lysis

Step 2.

A 4.5 mL volume of lysis buffer (20 mM Tris-HCl, pH 8.0; 100 mM NaCl; 50 mM Na₂EDTA), 20% SDS (250 μ L) and proteinase K (250 μ L, 10 mg/mL) were added to 1 g of frozen cells, which were then incubated for 2 h at 60 °C with intermittent mixing. Then 250 μ L of proteinase K (10 mg/mL) were added again, and the cells were incubated for 40 min under the same conditions.

Extraction

Step 3.

The lysate was mixed for 10 min with 8 mL of phenol saturated with TE buffer, pH 8.0 using the tube rotator. The extract was centrifuged (10 min, 16,100 g, 20 °C), then the aqueous phase was taken and again mixed with an equal volume of phenol: chloroform: isoamyl alcohol mixture (25:24:1). It was then mixed and centrifuged under the same conditions. The aqueous phase was mixed with a chloroform: isoamyl alcohol mixture (24:1), mixed and centrifuged under the same conditions. Aqueous phase was used for the subsequent stages.

DNA precipitation

Step 4.

DNA was precipitated by adding 1/10 volume of 3 M sodium acetate, pH 5.2, and 1 volume of isopropanol to the aqueous phase. The solution was matured at 0 °C for 1 h and then centrifuged (30 min, 16,100 g, 4

°C). The supernatant was removed, and the precipitate was twice rinsed with 70% ethanol, dried in air and dissolved in sterile TE buffer (10 mM Tris-HCl; 1 mM Na₂EDTA, pH 8.0).
