

# SYBR Gold working solutions

Li Deng

## Abstract

**Citation:** Li Deng SYBR Gold working solutions. **protocols.io**

dx.doi.org/10.17504/protocols.io.c7azid

**Published:** 04 Jan 2016

## Protocol

### Step 1.

Thaw the commercial stock of SYBR Gold in the dark at RT

### Step 2.

Centrifuge at 3000 rpm for 5 minutes

 **DURATION**

00:05:00

 **NOTES**

**VERVE Team** 24 Jun 2015

Because SYBR Gold is in DMSO.

### Step 3.

Dilute SYBR Gold in 0.02  $\mu\text{m}$  filtered TE buffer to 100x (10  $\mu\text{l}$  in 990  $\mu\text{l}$  TE buffer)

### Step 4.

0.02  $\mu\text{m}$  filter the diluted SYBR (100x)

 **NOTES**

**VERVE Team** 24 Jun 2015

This working stock can be stored at  $-20^{\circ}\text{C}$  and re-thawed one time.