

# Cvijl Purification From IL-3A Virus Infected NC64A Chlorella

David Dunigan and Irina Agarkova

## Abstract

**Citation:** David Dunigan and Irina Agarkova Cvijl Purification From IL-3A Virus Infected NC64A Chlorella. **protocols.io** dx.doi.org/10.17504/protocols.io.er3bd8n

**Published:** 13 Jun 2016

## Guidelines

### MATERIALS AND ASSAY CONDITIONS:

1) 7 hour IL-3A virus infected NC64A, pellets frozen at -80°C

2) Buffer A:

- 10 mM Tris-Acetate, pH 8.0
- 10 mM 2-M
- 50 µg/mL PMSF

3) Buffer B:

- 20 mM Tris-Acetate, pH 8.0
- 0.5 mM EDTA
- 7 mM 2-ME, 10% Glycerol

4) Buffer B, pH 8.5:

- 20 mM Tris-Acetate, pH 8.5
- 0.5 mM EDTA
- 7 mM 2-ME

5) Storage buffer:

- 20 mM Tris-Acetate, pH 8.0
- 0.5 mM EDTA
- 0.1 mM DTT
- 50 mM KOAc
- 5 mM MgAc, 50% Glycerol
- 100 µg/mL BSA (the BSA is added after dialysis into the storage buffer from a 10 mg/mL stock)

6) 4 M NaCl

7) 28% Polyethyleneglycol (PEG) 8000

8) Assay buffer and conditions:

1X Cvijl assay buffer:

- 20 mM GlycylGlycine, pH 8.5 (with KOH)
- 10 mM MgAc

- 0.1 mM DTT
- 50 mM KOAc
- 100 µg/mL casein (optional)

All assays are carried out in 20.0 µL volumes with 1 µg of pUC19 DNA as substrate for 60 to 120 min at 25°C. The assays are electrophoresed on 2.0% agarose gels (100 mL gels with double 20 lane combs) for 1 hour at 100 volts in 1X TPE buffer. Gels are stained with 0.5 µg/mL ethidium bromide for 30 min and photographed on a UV light box.

## Before start

7 hour IL-3A virus infected NC64A, pellets frozen at -80°C

## Protocol

### Step 1.

Prepare Buffer A:



[PROTOCOL](#)

#### . [Cvijl Buffer A](#)

CONTACT: [Irina Agarkova](#)

#### Step 1.1.

10 mM Tris-Acetate, pH 8.0

#### Step 1.2.

10 mM 2-ME

#### Step 1.3.

50 µg/mL PMSF

### Step 2.

Prepare Buffer B:



[PROTOCOL](#)

#### . [Cvijl Buffer B](#)

CONTACT: [Irina Agarkova](#)

#### Step 2.1.

20 mM Tris-Acetate, pH 8.0

#### Step 2.2.

0.5 mM EDTA

#### Step 2.3.

7 mM 2-ME, 10% Glycerol

### Step 3.

Prepare Buffer B, pH 8.5:



[PROTOCOL](#)

## . [Cvijl Buffer B, pH 8.5](#)

CONTACT: [Irina Agarkova](#)

### **Step 3.1.**

20 mM Tris-Acetate, pH 8.5

### **Step 3.2.**

0.5 mM EDTA

### **Step 3.3.**

7 mM 2-ME

### **Step 4.**

Prepare Storage Buffer:



## . [Cvijl Storage Buffer](#)

CONTACT: [Irina Agarkova](#)

### **Step 4.1.**

20 mM Tris-Acetate, pH 8.0

### **Step 4.2.**

0.5 mM EDTA

### **Step 4.3.**

0.1 mM DTT

### **Step 4.4.**

50 mM KOAc

### **Step 4.5.**

5 mM MgAc, 50% Glycerol

### **Step 4.6.**

100 µg/mL BSA (the BSA is added after dialysis into the storage buffer from a 10 mg/mL stock)

### **Step 5.**

Prepare 1X Cvijl Assay Buffer:



## . [1X Cvijl Assay Buffer](#)

CONTACT: [Irina Agarkova](#)

### **Step 5.1.**

20 mM GlycylGlycine, pH 8.5 (with KOH)

### **Step 5.2.**

10 mM MgAc

### **Step 5.3.**

0.1 mM DTT

### **Step 5.4.**

50 mM KOAc

### **Step 5.5.**

100 µg/mL casein (optional)

**Step 6.**

Thaw 7 hour IL-3A virus infected NC64A *chlorella* and suspend in MSK flasks with Buffer A.

**📌 NOTES**

**Irina Agarkova** 29 Mar 2016

Suspend with 20 mL per flask per  $1.0\text{--}1.5 \times 10^{11}$  infected cells.

**Step 7.**

Homogenize the cells in the MSK mechanical homogenizer with 15 gm of 0.3 mm glass beads at 4,000 rpm for 90 sec (2 X 45 sec) with CO<sub>2</sub> cooling.

**🕒 DURATION**

00:01:30

**Step 8.**

Recover the homogenate to clean tubes.

**Step 9.**

Wash the glass beads 3X with 5 mL of Buffer A and combine with the homogenate.

**Step 10.**

Centrifuge the homogenate in the Sorvall SS34 rotor at 10,000 rpm, 20 min, 4°C.

**🕒 DURATION**

00:20:00

**Step 11.**

Save the supernatant.

**Step 12.**

Adjust the homogenate supernatant to 70% saturation with (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> at 4°C with gentle stirring.

**📌 NOTES**

**Irina Agarkova** 29 Mar 2016

Add the (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> gradually.

**Step 13.**

Incubate at 4°C for 60-90 min without stirring.

**🕒 DURATION**

01:30:00

**Step 14.**

Centrifuge the material in the Sorvall SS34 rotor at 10,000 rpm 10 min, 4°C.

**🕒 DURATION**

00:10:00

**Step 15.**

Save the pellet.

**Step 16.**

Suspend the pellets with Buffer A.

**Step 17.**

Per mL of suspension add: 0.45 mL of 4 M NaCl and 0.45 mL of 28% PEG 8000 (heated to 65°C).

**Step 18.**

Mix gently by inversion for 5-10 min.

**🕒 DURATION**

00:10:00

**Step 19.**

Centrifuge the material in the Sorvall SS34 rotor at 10,000 rpm, 10 min, 4°C.

## DURATION

00:10:00

### Step 20.

Save the supernatant.

## NOTES

**Irina Agarkova** 29 Mar 2016

If there is small particulate matter in the supernatant, centrifuge a second time as before and save the supernatant.

### Step 21.

Dilute the supernatant with 10-15 volumes of Buffer B to reduce the NaCl concentration.

### Step 22.

Load the material overnight onto a Heparin-Sepharose column equilibrated with Buffer B in the cold room.

## DURATION

18:00:00

### Step 23.

Elute the Heparin-Sepharose column with Buffer B using a 0-2.0 M KOAc gradient.

### Step 24.

Assay the column fractions and pool the active fractions.

## PROTOCOL

### . [Cvijl Assay Conditions](#)

CONTACT: [Irina Agarkova](#)

#### Step 24.1.

All assays are carried out in 20.0  $\mu$ L volumes with 1  $\mu$ g of pUC19 DNA as substrate for 60 to 120 min at 25°C.

## DURATION

02:00:00

#### Step 24.2.

The assays are electrophoresed on 2.0% agarose gels (100 mL gels with double 20 lane combs) for 1 hour at 100 volts in 1X TPE buffer.

## DURATION

01:00:00

#### Step 24.3.

Gels are stained with 0.5  $\mu$ g/mL ethidium bromide for 30 min and photographed on a UV light box.

## DURATION

00:30:00

### Step 25.

Dilute the pooled fractions with 10-15 volumes of Buffer B to reduce the salt concentration.

### Step 26.

Load the material overnight onto a Blue-Sepharose column equilibrated with Buffer B in the cold room.

## DURATION

18:00:00

**Step 27.**

Elute the Blue-Sepharose column with Buffer B using a 0-2.0 M KOAc gradient.

**Step 28.**

Assay the column fractions and pool the active fractions.

**. [CviJI Assay Conditions](#)**

CONTACT: [Irina Agarkova](#)

**Step 28.1.**

All assays are carried out in 20.0 µL volumes with 1 µg of pUC19 DNA as substrate for 60 to 120 min at 25°C.



02:00:00

**Step 28.2.**

The assays are electrophoresed on 2.0% agarose gels (100 mL gels with double 20 lane combs) for 1 hour at 100 volts in 1X TPE buffer.



01:00:00

**Step 28.3.**

Gels are stained with 0.5 µg/mL ethidium bromide for 30 min and photographed on a UV light box.



00:30:00

**Step 29.**

Dilute the pooled fractions with 10-15 volumes of Buffer B, pH 8.5 to reduce the salt concentration.

**Step 30.**

Load the material overnight onto a Q-Sepharose column equilibrated with Buffer B, pH 8.5 in the cold room.



18:00:00

**Step 31.**

Elute the Q-Sepharose column with Buffer B, pH 8.5 using a 0-2.0 M KOAc gradient.

**Step 32.**

Assay the column fractions and pool the active fractions.

**. [CviJI Assay Conditions](#)**

CONTACT: [Irina Agarkova](#)

**Step 32.1.**

All assays are carried out in 20.0 µL volumes with 1 µg of pUC19 DNA as substrate for 60 to 120 min at 25°C.



02:00:00

**Step 32.2.**

The assays are electrophoresed on 2.0% agarose gels (100 mL gels with double 20 lane combs) for

1 hour at 100 volts in 1X TPE buffer.

🕒 DURATION

01:00:00

**Step 32.3.**

Gels are stained with 0.5 µg/mL ethidium bromide for 30 min and photographed on a UV light box.

🕒 DURATION

00:30:00

**Step 33.**

Concentrate the pooled enzyme by dialysis overnight into storage buffer at 4°C.

🕒 DURATION

18:00:00

**Step 34.**

Add BSA (10 mg/mL) to a final concentration of 100 µg/mL.

**Step 35.**

Store the enzyme at -20°C.