



Nov 15,  
2019

## BDA Histology

Blake Butler<sup>1</sup>

<sup>1</sup>University of Western Ontario

1

Works for me

dx.doi.org/10.17504/protocols.io.z4af8se



Blake Butler  
University of Western Ontario



### ABSTRACT

Staining protocol for visualization of BDA 3K.

Species - Feline

Tissue Preparation - Perfusion of 4% PAF; Cryoprotected 30% Sucrose; Section thickness of 60 microns

### STEPS MATERIALS

NAME	CATALOG #	VENDOR
Avidin/Biotin HRP Complex	PK-6100	Vector Laboratories
Normal Goat Serum	CL1200-500	Cedarlane
3,3'-Diaminobenzidine tetrahydrochloride	D5905	Sigma Aldrich

### Preparatory Steps

- 1 Reconstitute the Vectastain ABC solution

Prepare 80 mL solution:

79 ml [M] 0.1 Molarity (M) Physiological Buffer

480 µl Solution A

480 µl Solution B



Avidin/Biotin HRP Complex

by Vector Laboratories

Catalog #: PK-6100

Thaw 6 x 1mL vials of NGS for later use



Normal Goat Serum

by Cedarlane

Catalog #: CL1200-500

## Histological Reaction

2

2.1 Rinse sections in 0.1 M Physiological Buffer (4x5min)

🕒 00:05:00

🕒 00:05:00

🕒 00:05:00

🕒 00:05:00

2.2 Block Endogenous Peroxidase

🧴 105 ml Methanol

🧴 45 ml De-ionized Water

🧴 750 µl (M) 30 % volume Hydrogen Peroxide

🕒 00:30:00

2.3 Rinse sections in 0.1 M Physiological Buffer (4x5min)

🕒 00:05:00

🕒 00:05:00

🕒 00:05:00

🕒 00:05:00

2.4 Transfer tissue to individual reaction wells and incubate in Avidin/Biotin HRP Complex Prepared in Step 1

🕒 01:30:00

2.5 Transfer tissue back to reaction trays and rinse sections in 0.1 M Physiological Buffer (3x10min)

🕒 00:10:00

🕒 00:10:00

🕒 00:10:00

## 2.6 Incubate in DAB-Nickel Chromogen Solution

🧴 **150 ml** [M] **0.1 Molarity (M)** Physiological Buffer

🧴 **0.1125 g** 3,3' Diaminobenzidine Tetrachloride (DAB)

🧴 **1.875 ml** [M] **1 Mass / % volume** Cobalt Chloride (Added dropwise)

🧴 **3 ml** [M] **1 Mass / % volume** Nickel Ammonium Sulfate (Added dropwise)

🧴 **50 µl** [M] **30 % volume** Hydrogen Peroxide (Added 60s before incubation)



3,3'-Diaminobenzidine  
tetrahydrochloride

by Sigma Aldrich

Catalog #: D5905

🕒 **00:20:00** (Approximate time - Actual endpoint determined visually)

## 2.7 Rinse sections in 0.1 M Physiological Buffer (3x5min)

🕒 **00:05:00**

🕒 **00:05:00**

🕒 **00:05:00**

## 3 Mount from 0.01 M Physiological Buffer



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited