



Jan 17,
2020

High-Throughput Beta-glucuronidase (GUS) assay for *Phaeodactylum tricornutum*

Erin Garza¹, Vincent Bielinski²

¹J. Craig Venter Institute, ²J. Craig Venter Institute, Synthetic Biology & Bioenergy Group

1 Works for me [dx.doi.org/10.17504/protocols.io.bbexijfn](https://doi.org/10.17504/protocols.io.bbexijfn)

Protist Research to Optimize Tools in Genetics (PROT-G)

JCVI West Protocols

Erin Garza
J. Craig Venter Institute

ABSTRACT

A high-throughput method for measuring β -glucuronidase (GUS) activity in the diatom *Phaeodactylum tricornutum*. This protocol has been optimized for 250 μ l volumes. For larger volumes see the following protocol [dx.doi.org/10.17504/protocols.io.hefb3bn](https://doi.org/10.17504/protocols.io.hefb3bn), which this protocol was based off of.

MATERIALS

NAME	CATALOG #	VENDOR
MUG	MUG	Gold Biotechnology
Sodium carbonate	222321	Sigma – Aldrich
B-PER™ Bacterial Protein Extraction Reagent	78243	Thermo Fisher

MATERIALS TEXT

Flat bottom transparent and opaque 96-well plates

GUS extraction buffer- 50 mM NaPO₄H₂ (pH 7), 0.1% Triton X-100, and 10 μ M β ME + 1 mM 4-Methylumbelliferyl β -D-Glucuronide (MUG)

GUS stop buffer- 0.2 M Na₂CO₃

Plate reader

Swing bucket centrifuge with plate adapter

BEFORE STARTING

Phaeodactylum tricornutum cultures were initially grown in 5 ml L1 + antibiotics in a 50-ml conical at 18°C until the cell concentration reached at least 1x10⁶ cells ml⁻¹.

Centrifuge

- 1 Transfer 250 μ l of each *P. tricornutum* culture to a 96-well plate and centrifuge at 3000 x g for 10 min. Discard supernatant.

Lyse

- 2 To lyse the cells, add 150 μ l bacterial protein extraction reagent (B-PER, ThermoFisher) to each well and mix by pipetting.

Centrifuge

- 3 Centrifuge plate for 10 min at 3000 x g. Transfer supernatants to a new 96-well plate, being careful not to disturb the cell debris.

Extract

- 4 Transfer 50 µl of each lysate to a new plate and add 125 µl GUS extraction buffer + 1 mM MUG to each well. Incubate the plate for 1 h at 37°C.

GUS extraction buffer= 50 mM NaPO₄H₂ (pH 7), 0.1% Triton X-100, and 10 µM βME

Stop Reaction

- 5 To stop the reaction, add 150 µl GUS stop buffer (0.2 M Na₂CO₃) to each well and mix by pipetting.
- 6 Transfer 200 µl quenched reaction to an opaque 96-well plate.

Read fluorescence

- 7 Determine fluorescence using a plate reader. Settings: excitation- 360 nm; emission- 440 nm.
- 8 If fluorescence readings are too high to get a readout, dilute with additional stop buffer.

Normalization

- 9 Use remaining cell lysates to perform a BCA assay (or an equivalent assay) to normalize the GUS activity to total cell protein for each culture.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited