

Oct 20, 2019

## Cell-ELONA

Guillermo Fernández Rodríguez<sup>1</sup><sup>1</sup>AEGIS - Madrid iGEM 2019
1 Works for me [dx.doi.org/10.17504/protocols.io.7umhnu6](https://doi.org/10.17504/protocols.io.7umhnu6)

AEGIS - Madrid iGEM 2019

Guillermo Fernández Rodríguez ⚡

## ABSTRACT

ELONA (Enzyme-Linked Oligonucleotide Assay), is a biochemical method based on enzyme linked immunosorbent assay (ELISA) that allows to demonstrate that aptamers selected by SELEX can be effective and useful as biorecognition molecules and laboratory tools. The ELONA format choosed uses an anti-digoxigenin antibody to recognize an aptamer previously labelled with digoxigenin. This antibody is conjugated with a peroxidase enzyme, and once it adds ABTS solution, it will be responsible for the colourimetric reaction which will be detected.

(We used 5 replicates per dilution)

## MATERIALS



NAME ▾	CATALOG # ▾	VENDOR ▾
MG	26528	addgene
Anti-Digoxigenin-AP, Fab fragments	11093274910	Sigma – Aldrich
Sodium bicarbonate	S6014	Sigma – Aldrich
PBST (PBS 1:1000 Tween-20)	<a href="#">View</a>	
BSA		Sigma Aldrich
Centrifuge	5415D	Eppendorf Centrifuge
LB	L24400-2000.0	Research Products International (rpi)
ABTS, solution	002024	Thermo Fisher
White 96-Well Immuno Plates, Maxisorp, Flat-Bottom, MaxiSorp, 350µL	436110	Thermo Fisher

## BEFORE STARTING



Clean all the working surface with ethanol.

## Pre-Coating











- 1 Inoculate a single colony of E.Coli DH5α from LB agar plate in 10 ml of LB. Use a sterile pipette tip, selecting a single colony from LB agar plate. The liquid culture is incubated overnight 37 °C
- 2 Spin at 4000 rpm for 00:05:00 . Discard the supernatant, collect pellet and re-suspend in 10 ml of NaHCO<sub>3</sub>-Na<sub>2</sub>CO<sub>3</sub>, 50mM, pH 9,6. Mix by inverting the tube.

- 3 Spin at 4000 rpm for  00:05:00 . Discard the supernatant, collect pellet and re-suspend in  8 ml of NaHCO<sub>3</sub>-Na<sub>2</sub>CO<sub>3</sub>, 50mM, pH 9,6. Mix by inverting the tube.
- 4 Read the absorbance (600nm). Dilute the sample with NaHCO<sub>3</sub>-Na<sub>2</sub>CO<sub>3</sub>, 50mM, pH 9,6, and adjust the absorbance to 1.

#### Coating (Automated)

- 5 Dilute the sample 1:4 with NaHCO<sub>3</sub>-Na<sub>2</sub>CO<sub>3</sub>, 50mM, pH 9,6 in a eppendorf tube.
- 6 Add  200 µl of the bacterial sample into Nunc MaxiSorp 96-well plate and Coating buffer for the negative control. Incubate overnight at  4 °C with agitation (260rpm).

#### Cell-ELONA (Automated)

- 7 Remove Coating buffer. CAUTION: Be careful to do not touch the well and remove the bacteria.
- 8 Wash 3x  200 µl with PBS 1x-Tween 0,1%. Remove the drops after the last wash. CAUTION: Be careful to do not touch the well and remove the bacteria.
- 9 Block the plate with  200 µl PBS 1x BSA 5% for  01:00:00 (260 rpm).
- 10 Structure digoxigenin-labelled aptamers (denatured for 10 min at 95°C and then cooled for 10 min on ice), previously prepared in different concentrations.
- 11 Wash 3x  200 µl with PBS 1x tween 0,1%. Remove the drops after the last wash. CAUTION: Be careful to do not touch the well and remove the bacteria.
- 12 Add  100 µl /well of the structured aptamers. Incubated for  01:00:00 .
- 13 Wash 3x  200 µl with PBS 1x tween 0,1%. Remove the drops after the last wash.
- 14 Add anti-body antidigoxigenin (  100 µl /well) preparing 1:1000 dilution in Aptamer buffer-BSA 0,2%. Incubate at room temperature for  01:00:00 .
- 15 Wash 3 x  200 µl with PBS 1x Tween 0,1%.

- 16 Add 100  $\mu$ L/well of ABTS solution. Read the absorbance (405 nm) every 🕒 00:10:00 for 🕒 01:00:00 . ADVICE: We recomend you to buy the ABTS than comes diluted and with the oxygene peroxide.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited