



Oct 18, 2019

Promega BacTiter-Glo Assay

NUS iGEM¹

¹National University of Singapore

1 Works for me dx.doi.org/10.17504/protocols.io.8ekhtcw



NUS iGEM
National University of Singapore

MATERIALS

NAME

CATALOG #

VENDOR

BacTiter-Glo Microbial Cell Viability Assay

G8230

Promega

- 1 Prepare BacTiter-Glo reagent mix by transferring BacTiter-Glo buffer into BacTiter-Glo substrate



These reagents are part of a set - Promega's BacTiter-Glo Assay.

- 2 Aliquot 50 μ l of cells into two wells of 96-well plate (duplicates)
- 3 Add 50 μ l of BacTiter Glo reagent mix into each well
- 4 Include PBS-only wells as blank
- 5 Measure luminescence using microplate reader



Cytation 5

Automated microscope; Multi-mode plate reader

BioTek N.A. [Link](#)



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited