

# Complete Yeast Protein Extraction

Tobias von der Haar

## Abstract

This protocol is suitable for rapidly generating extracts of fungal cells for gel electrophoresis or western blots. It was optimised for complete extraction of difficult-to-solubilise proteins. This protocol was originally described in [von der Haar 2007](#).

**Citation:** Tobias von der Haar Complete Yeast Protein Extraction. **protocols.io**

dx.doi.org/10.17504/protocols.io.g8wbzxe

**Published:** 05 Mar 2017

## Protocol

### Step 1.

Prepare the basic solutions:

**Extraction Buffer:** 0.1 M NaOH, 0.05 M EDTA, 2% SDS

(Note: before use, this buffer is supplemented with 2% 2-mercaptoethanol. With 2-ME the shelf life is several hours, without 2-ME the shelf life of this solution is 1 week).

**Caution:** *This buffer is strongly alkaline and additional protective equipment (gloves and goggles) should be worn while preparing and handling it.*

**Neutralisation Buffer:** 4 M acetic acid in water. Neat acetic acid is 11 M in concentration.

**Loading buffer:** 0.25 M Tris-HCl pH 6.8, 50% glycerol, 0.05% bromophenolblue.

### Step 2.

Pre-heat a heating block to 96°C. The extraction works efficiently for most proteins at temperatures of 80°C or above.

### Step 3.

Harvest  $10^8$  yeast cells (about 5  $\text{oD}_{600}$  units of logarithmically growing yeast) by centrifuging in a microcentrifuge at max speed for 1 minute. Completely remove the supernatant.

#### **Step 4.**

Freshly supplement 1 ml of extraction buffer with 20 µl 2-mercaptoethanol (scale this up for more than 5 extractions). Resuspend the pelleted yeast cells in 200 µl of the supplemented extraction buffer. Incubate the resuspended cells at 96°C for 10 minutes.

#### **Step 5.**

Add 5 µl of neutralisation buffer, vortex the cells.

#### **Step 6.**

Add 50 µl of loading buffer, vortex again.

#### **Step 7.**

The extract is now ready to be loaded onto SDS-PAGE gels. Clearing by centrifugation can be done but is not required, and loading the uncleared (but well resuspended) cell slurry into the pockets of the gel aids extraction.

### **Warnings**

The extraction buffer is strongly alkaline. Additional protective equipment (gloves and goggles) should be worn while handling this buffer.