

# Proteomics Analysis

Audrey Dubourg, Dong Xia, John P Winpenny, Suha Al Naimi, Maha Bouzid, Darren W Sexton, Jonathan M Wastling, Paul R Hunter, Kevin M Tyler

## Abstract

This protocol provides an efficient and standardized way to prepare samples ready for mass spec analysis.

**Citation:** Audrey Dubourg, Dong Xia, John P Winpenny, Suha Al Naimi, Maha Bouzid, Darren W Sexton, Jonathan M Wastling, Paul R Hunter, Kevin M Tyler Proteomics Analysis. **protocols.io**

dx.doi.org/10.17504/protocols.io.mcec2te

**Published:** 04 Jan 2018

## Before start

Prepare solutions as described in the steps.

## Protocol

### Sample preparation

#### Step 1.

Dilute the amount of protein needed, taking into account complexity, into 25mM AmBic. Small Eppendorf tube preferred.

### Detergent treatment

#### Step 2.

Add 10µL of 1% (w/v) RapiGest (0.05% (w/v) final).

#### AMOUNT

10 µl Additional info: 1% (w/v) RapiGest

### Detergent treatment

#### Step 3.

Heat at 80°C, 10minutes, vortex briefly at 5min.

#### TEMPERATURE

80 °C Additional info: Heat

### Detergent treatment

#### Step 4.

Spin quickly to return liquid to the bottom of the tube.

#### Reduction

##### Step 5.

Add 10 µL of a 9.2mg/mL solution of DTT (3 mM final). Vortex mix.

#### AMOUNT

10 µL Additional info: 9.2mg/mL solution of DTT

#### Reduction

##### Step 6.

Incubate for 60°C, 10minutes.

#### TEMPERATURE

60 °C Additional info: Incubation

#### Reduction

##### Step 7.

Cool to RT and quickly spin to return liquid to the bottom of the tube.

#### Alkylation

##### Step 8.

Add 10 µL of a 33mg/mL solution of iodoacetamide (9 mM final). Vortex.

#### AMOUNT

10 µL Additional info: 33mg/mL solution of iodoacetamide

#### Alkylation

##### Step 9.

Incubate at RT, **IN THE DARK** for 30min.

#### Digestion

##### Step 10.

Add trypsin to 50:1 protein:trypsin ratio. Incubate 12-16h (overnight) at 37°C.

#### TEMPERATURE

37 °C Additional info: Incubation

#### Mass spec analysis

##### Step 11.

Analyze peptide mixtures on either nanoACQUITY-nLC system (Waters MS technologies) or Ultimate 3000 nano system (Thermo Fisher Scientific) followed by an LTQ-Orbitrap Velos (ThermoFisher Scientific) mass spectrometer or a Q-Exactive mass spectrometer (Thermo Fisher Scientific).