

# Chemically competent *E. coli* cells

Carlos Helbig

## Abstract

This protocol describes how to make chemically competent *Escherichia coli* cell.

For the protocol *E. coli* strain NEB Turbo was taken.

**Citation:** Carlos Helbig Chemically competent *E. coli* cells. **protocols.io**

[dx.doi.org/10.17504/protocols.io.psydnfw](https://dx.doi.org/10.17504/protocols.io.psydnfw)

**Published:** 30 Apr 2018

## Before start

Make sure you have the following solutions:

- **TYM**

- Bacto-tryptone 20g/L
- Yeast extract 5g/L
- NaCl 100mM
- MgSO<sub>4</sub> 10mM

- **TfB1**

- KAc 30mM (pH5.8)
- KCl 100mM
- CaCl<sub>2</sub> 10mM
- Glycerol 15%
- (Add after autoclaving 5% vol. MnCl<sub>2</sub> 1M)

- **TfB2**

- Na-MOPS 10mM
- CaCl<sub>2</sub> 75mM
- KCl 10mM
- Glycerol 15%
- Adjust pH 6.8 - 7.0

## Protocol

### Step 1.

Grow 50ml over night culture in TYM medium

**Step 2.**

Transfer 10ml cells to 500ml TYM medium

**Step 3.**

Incubate: OD=0.55, 37°C, 200rpm

**Step 4.**

Cool cells on ice

**Step 5.**

Centrifuge: 3500g, 15min., 4°C

**Step 6.**

Remove supernatant completely

**Step 7.**

Resuspend pellets in 100ml ice cold Tfb1 buffer

**Step 8.**

Centrifuge: 3500g, 15min., 4°C

**Step 9.**

Resuspend pellets in 30ml ice cold Tfb2 buffer

**Step 10.**

Aliquot the cells into chilled tubes (50µl aliquots)

**Step 11.**

Freeze the cells in liquid nitrogen

**Step 12.**

Store at -80°C