



Oct 06, 2019

## Haematoxylin and eosin staining

Norfitriah Mohamed Sohaimi<sup>1</sup>, Mohd Hair Bejo<sup>1</sup>, Abdul Rahman Omar<sup>1</sup>, Aini Ideris<sup>1</sup>, Nurulfiza Mat Isa<sup>2</sup>

<sup>1</sup>Faculty of Veterinary Medicine, Universiti Putra Malaysia, Serdang, Selangor, Malaysia., <sup>2</sup>Faculty of Biotechnology and Biomolecular Sciences, Universiti Putra Malaysia, Serdang, Selangor, Malaysia.

1 Works for me dx.doi.org/10.17504/protocols.io.7zghp3w

Norfitriah Mohamed Sohaimi

### MATERIALS

| NAME ▾                                    | CATALOG # ▾ | VENDOR ▾      |
|---|-------------|---------------|
| Dulbecco's Modified Eagle's Medium (DMEM) | D5796       | Sigma Aldrich |

### MATERIALS TEXT

For detection of intranuclear inclusion bodies, primary CEL cells were grown on sterile coverslips [Adair, 1978]. Briefly, the growth medium was aspirated out from confluent monolayer and washed twice with sterile medium free serum prior to inoculation of 0.1mL viral supernatant. Culture was added with maintenance medium, DMEM and 2% FBS and kept for 48 hours incubation, harvested and fixed in 10% buffered formalin for 5 minutes and stained with haematoxylin and eosin (HE).



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited