



Feb 27, 2019

Working

iPSC gDNA Extraction: For Screening Edited Clones

In 1 collection

Celeste Karch¹, Rita Martinez¹, Jacob Marsh¹¹Washington University in St Louis

dx.doi.org/10.17504/protocols.io.x8afrse

Neurodegeneration Method Development Community

Tech. support email: ndcn-help@chanzuckerberg.com

Celeste Karch

Washington University in St Louis



Comprehensive Genomic
Editing and Screening
Protocol Updated
02142019.docx

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

This protocol is part of the [Screening Edited iPSC Clones collection](#).

SAFETY WARNINGS

Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.

1 Spin 96 well PCR plate at 3800 rpm for 00:30:00 at 4 °C .

2 Pipet off supernatant.



If worried about removing cells, transfer supernatant with multichannel pipet into new PCR plates and store at -20 °C until DNA extraction is complete.

3 To cell pellet (often not visible), add 50 µl QuickExtract DNA Solution ([Epicentre Technologies QE09050](#)).

4 Vortex plate for 00:00:15 .

- 5 Incubate plate at  **65 °C** for  **00:06:00** .
- 6 Vortex plate for  **00:00:15** .
- 7 Incubate plate in thermocycler at  **98 °C** for  **00:02:00** .
- 8 Place plate in  **-20 °C** for storage, until ready to use for further screening experiments .



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited