



Dec 09, 2019

SPRI purification beads

Tomasz Suchan¹¹W. Szafer Institute of Botany, Polish Academy of Sciences

1

Works for me

dx.doi.org/10.17504/protocols.io.q3qdywmw

Molecular Biogeography Group



Tomasz Suchan

W. Szafer Institute of Botany, Polish Academy of Sciences



ABSTRACT

Protocol for preparing SPRI beads 100x cheaper than the commercial AMPure XP.

Based on DeAngelis et al. (1995) and online protocols:

https://s3-us-west-2.amazonaws.com/oww-files-public/1/17/Serapure_bead_recipe.pdfhttps://s3-us-west-2.amazonaws.com/oww-files-public/f/f8/SPRI_buffers_v2_2.pdf

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

DeAngelis, M. M., Wang, D. G., & Hawkins, T. L. (1995). Solid-phase reversible immobilization for the isolation of PCR products. *Nucleic acids research*, 23(22), 4742.

MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

PEG-8000

NaCl

Tween-20

P-7949

Sigma-aldrich

Water, nuclease free

1 M Tris-HCl pH 8.0

0.5 M EDTA pH 8.0

Sera-Mag Carboxylate-Modified Magnetic Particles (Hydrophylic), 15 mL

24152105050250

Ge Healthcare

Preparation of the solutions

- 1 Prepare 2.5 M solution of NaCl:
 - place a 50 ml Falcon tube on the scale and tare
 - add 14.61 g of NaCl to the tube
 - add water to 40 ml
 - dissolve NaCl
 - fill up the tube with water to 50 ml and mix

- 2 Prepare 50% solution of PEG-8000:
 - place a 50 ml Falcon tube on the scale and tare
 - add 12.5 g of PEG-8000 to the tube
 - add 14 ml of nuclease-free water
 - close the tube, and dissolve PEG by mixing vigorously, preferably on an orbital shaker with heating
 - fill up the tube with water to 25 ml and mix gently by flipping the tube to avoid bubbles
- 3 Prepare 10% Tween-20:
 1. add 900 ul of water to 1.5 ml tube
 2. slowly aspirate 100 ul of Tween-20 with just a tip of the pipette submerged to avoid liquid retention outside of the tip (it is very viscous so take time to aspirate the whole volume)
 3. add Tween-20 to the water, mix by pipetting
- 4 Prepare TE buffer by mixing:
 - 49.4 ml of molecular-grade water
 - 500 ul of 1M Tris solution
 - 100 ul of 0.5M EDTA solution



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited