



May 14,
2019

Working

UC Davis - IN-VIVO Glucose-stimulates Insulin Secretion Test [↗](#)

Fawaz G. Haj¹

¹University of California, Davis

[dx.doi.org/10.17504/protocols.io.yptfvnn](https://doi.org/10.17504/protocols.io.yptfvnn)

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

An in-vivo glucose-stimulates insulin secretion test is designated to determine alterations in insulin secretion by the pancreas upon a bolus IP Glucose injection.

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=86&docType=Protocol>

MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾	CAS NUMBER ▾	RRID ▾
45% Glucose solution	NC0025179(50-165-7017 replaced)	Fisher Scientific		
Insulin Syringes	14-826-79	Fisher Scientific		
Saline Solution	L97753	Fisher Scientific		
Ultra Sensitive Mouse Insulin ELISA kit	90080	Crystal chem		AB_2783626
Heparin Sodium	401586B	Abraxis		
Easy Check Glucose test strips	00-101(new SKU 88982400)	JRS Medical		
Easy Check Glucose monitor	Y4209 (new SKU 88972401)	JRS Medical		

MATERIALS TEXT

Dilute the glucose stock solution (45%) with saline to 20% by adding 20ml stock to 25ml 0.9% (w/v) sterile saline.

Note:

Fisher Scientific, [RRID:SCR_008452](#)

Ultra Sensitive Mouse Insulin ELISA kit, Cite this, ([Crystal Chem Cat# 90080](#), [RRID:AB_2783626](#))

- 1 Fast mice for 16 hours by taking away food the day before (3:00pm)
- 2 The following day, Calibrate the glucose meter according to the manufacturer's instructions.

Deprive mice from water then measure blood glucose level using a glucometer and remove immediately approximately 50 µl of blood from

- 3 the tail via a tail tip cut and transfer directly onto a sterile 0.5ml microcentrifuge tube containing 2ul Heparin.
- 4 Centrifuge at 8000rpm for 5min then transfer plasma (supernatant) to a new 0.5ml microcentrifuge tube and freeze at -80°C.
- 5 Give the mouse an intraperitoneal injection of Glucose (2g/kg) with a 27 G needle
- 6 Continue to take blood samples from the initial tail cut at 2, 5, 15 and 30 min flowing injection and repeat step3 and 4.

7 **NOTE:**

At the end of the experiment, wipe tail with 70% alcohol and allow drying. Ensure that blood loss from the tail stopped before placing the animal back to its cage.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited