



Oct 14, 2019

## Sanger sequencing

[iGEM Dusseldorf<sup>1</sup>](#)<sup>1</sup>Heinrich-Heine Universität Düsseldorf

1

Works for me

[dx.doi.org/10.17504/protocols.io.76phrdn](https://doi.org/10.17504/protocols.io.76phrdn)

iGEM Dusseldorf ⚡

### ABSTRACT

1. Label 1.5 mL tubes with barcodes from GATC Biotech (stored on top shelf in our office)
2. Add template (400 to 500 ng), primer and water (10 uL total volume, details below)

	uL per Reaction		Barcode	Template	Primer
Template	#ERROR!				
10 uM Primer	2.5				
Water	#ERROR!				
Total Volume	10.0				

3. Place tubes inside a small plastic bag and drop off in GATC Collection Box by 5:30 pm.
4. Results are usually ready to download from the [gatc-biotech.com](https://gatc-biotech.com) or Dropbox (details below) the next business day after 12 pm. The GATC collection point is located at the exit of building 26.43, opposite to building 26.53. Daily pick up of your samples from Monday to Friday 5:30 pm.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited