

# Diet-induced obesity murine model

Bianca Hemmeryckx, Dries Bauters, H. Roger Lijnen

## Abstract

Obesity can be induced by exposing mice to a western-type diet high in sugar and fat for 15 weeks. Control animals receive a standard fat diet for 15 weeks.

**Citation:** Bianca Hemmeryckx, Dries Bauters, H. Roger Lijnen Diet-induced obesity murine model. **protocols.io**

dx.doi.org/10.17504/protocols.io.kbacsie

**Published:** 15 Oct 2017

## Guidelines

-HFD is stored at +4°C

-Mouse needs to be out of cage, relaxed on a bench, a scalpel is used to cut off 1 mm of tail, the tail starts to bleed, glucose levels are measured on the second drop of blood.

## Protocol

### Step 1.

Randomize 8-weeks old male mice (C57BL/6J, C57BL/6N or gene-deficient mice with control littermates) to two groups based on body weight.

### Step 2.

Give one group 7g standard chow (SFD, 10.9 kJ/g) each day/mouse, and the other group 7g western high-fat diet (HFD, 22 kJ/g; kcal from fat 42%, 43% from carbohydrates and 15% from protein) for 15 weeks. As food intake averages 4g/day, the animals are fed ad libitum. Provide drinking water ad libitum.



### REAGENTS

 ssniff® EF R/M acc. TD88137 mod E15721-34 by Contributed by users

### Step 3.

House the animals in a temperature (20-23°C) and humidity (30-40%)-controlled room with a 12-hour light/dark cycle (7AM/7PM). Give 1-2 small wooden sticks to every cage containing mice on a HFD in order to allow correct wearing off of their teeth.

### Step 4.

Replace drinking water and change the cages every week. Check the health of the animals every day.

**Step 5.**

Replace each diet weekly, record food intake and body weight weekly.

**Step 6.**

After 15 weeks of diet administration, fast animals for 6 hours (7:30 AM – 1:30 PM).

**Step 7.**

Determine blood glucose levels on a drop of blood taken from the nicked tail of a free-moving mouse using the Accu-Chek Aviva meter and Accu-Chek Aviva glucose strips.

**REAGENTS**

- ✓ Accu-Chek Aviva glucose strips 06453970 by Contributed by users

**Step 8.**

Sacrifice animals by intraperitoneal injection of 60 mg/kg sodium pentobarbital (Nembutal) or Dolethal.

-Nembutal (60 mg/ml): we make 10x dilution with saline.

-Dolethal (200 mg/ml): we make 30x dilution with saline.

**REAGENTS**

- ✓ Nembutal 60 mg/ml by Contributed by users
- ✓ Dolethal 200 mg/ml by Contributed by users

**Step 9.**

When animals are asleep, collect blood from the retro-orbital sinus on trisodium citrate and processed for plasma (<http://dx.doi.org/10.17504/protocols.io.j9dcr26>). Store plasma in aliquots at -80°C.

**Step 10.**

Immobilize the animals, cut open the chest and wash out the heart and the entire circulatory system with saline by putting a perfusion needle in the apex (left ventricle). Make a small cut in the right atrium to allow blood to flow out.

**REAGENTS**

- ✓ Microflex Infusion set, orange 240.05 by Contributed by users

**Step 11.**

Make a midline incision in the skin and make a incision at the height of the hind limbs.

**Step 12.**

Pull skin to the side so the inguinal subcutaneous white adipose tissue (sWAT) depots are visible. Remove the internal lymph nodes. Dissect out the sWAT depots.

**Step 13.**

Make a midline incision in the muscle layer to open the abdominal cavity.

**Step 14.**

Locate the intra-abdominal gonadal (GON) WAT depots and liver. Removed them.

**Step 15.**

Flip the mouse, locate the interscapular brown adipose tissue depot and remove it.

**Step 16.**

Weigh all dissected organs.

**Step 17.**

Make portions of each tissue for RNA, protein extraction or histology.

**Step 18.**

Snap-freeze the portions for RNA and protein extraction and store at -80°C.

**Step 19.**

Process the portions for histology (see separate protocol).

**Step 20.**

Clean dissection material and bench with Vircon S, detergent, water and 70% ethanol.

**REAGENTS**

Virkon™ S 50 tablets/bottle NC9821357 by [Fisher Scientific](#)

**Step 21.**

Dispose of the animal body according to internal institutional guidelines.