

Zymo Research DNA Clean & Concentrator

Alan Cone

Abstract

Cleans and concentrates any DNA. Typically useful after PCR, restriction digest, or anything where DNA was altered and now you need your newly altered DNA in a pure form.

Citation: Alan Cone Zymo Research DNA Clean & Concentrator. **protocols.io**

dx.doi.org/10.17504/protocols.io.egbbbsn

Published: 22 Jan 2016

Protocol

Step 1.

In a 1.5 ml microcentrifuge tube, add 2-7 volumes of DNA Binding Buffer to each volume of DNA sample (see table below). Mix briefly by vortexing.

Application	DNA Binding Buffer : Sample	Example
Plasmid, genomic DNA (>2 kb)	2 : 1	200 µl : 100 µl
PCR product, DNA fragment	5 : 1	500 µl : 100 µl
ssDNA ¹ (e.g. cDNA, M13 phage)	7 : 1	700 µl : 100 µl



REAGENTS

Zymo DNA Binding Buffer [D4003-1-25](#) by [Zymo Research](#)

Step 2.

Transfer mixture to a provided Zymo-Spin™ Column in a Collection Tube.

Step 3.

Centrifuge for 30 seconds. Discard the flow-through.

Step 4.

Add 200 µl DNA Wash Buffer to the column. Centrifuge for 30 seconds. Repeat the wash step.



AMOUNT

200 µl Additional info:



REAGENTS

Zymo DNA Wash Buffer [D4003-2-6](#) by [Zymo Research](#)

Step 5.

Repeat the wash step.



AMOUNT

200 µl Additional info:



REAGENTS

Zymo DNA Wash Buffer [D4003-2-6](#) by [Zymo Research](#)

Step 6.

Add ≥ 6 μ l DNA Elution Buffer or water directly to the column matrix and incubate at room temperature for one minute. Transfer the column to a 1.5 ml microcentrifuge tube and centrifuge for 30 seconds to elute the DNA.

AMOUNT

6 μ l Additional info:

REAGENTS

Zymo DNA Elution Buffer [D3004-4-1](#) by [Zymo Research](#)

NOTES

Alan Cone 22 Jan 2016

After starting with 1 ug DNA, typical yields are between 50-100 ng / uL

Great recovery is 80%, typical recovery may be closer to 60%.