



Sep 20, 2019

High Resolution Imaging Mass Spectrometry Analysis using Bruker Daltonics Platforms

Elizabeth Neumann¹, Jamie Allen¹, David Anderson¹, Danielle Gutierrez¹, Jeff Spraggins¹

¹Vanderbilt University

1 Works for me dx.doi.org/10.17504/protocols.io.7gdhjs6

VU Biomolecular Multimodal Imaging Center

Human BioMolecular Atlas Program (HuBMAP) Method Development Community



Jamie Allen ⚡

ABSTRACT

Scope:

Acquire 10 μm spatial resolution imaging mass spectrometry (IMS) datasets of $\sim 5\text{mm}^2$ regions of tissue.

Expected Outcome:

A series of ion images that can be used to visualize the distribution of lipids localizing to physiological regions within the human kidney.

MATERIALS TEXT

Slide Adapter for Instrument

Bruker MALDI timsTOF Flex

BEFORE STARTING

Prepare tissue by following protocols for sectioning and matrix application.

- 1 Scann a 3200 dpi image of the matrix coated tissue section in the MTP slide adapter 2 with sufficient contrast to visualize the tissue boundaries.
- 2 Choose a method that has "height correct" initiated upon sample docking.
- 3 Insert the slide adapter plate into the instrument airlock. Then press "load."
- 4 While the height correct profile is being generated, open FlexImaging and follow the prompts to designate sample file name, method, pitch, scanned image of the slide, etc.
- 5 Once the height correct profile has been generated, train the target position with three teaching points following the instructions within FlexImaging.
- 6 Draw a measurement region around the tissue for IMS analysis.
- 7 Go to an area off of the tissue and ensure the focus of the laser is adequate for 10 μm spatial resolution imaging.

8 Start the acquisition.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited