

TAE buffer (1L of 50X)

Bao Thai

Abstract

Citation: Bao Thai TAE buffer (1L of 50X). **protocols.io**

dx.doi.org/10.17504/protocols.io.mpcc5iw

Published: 16 Jan 2018

Protocol

Step 1.

Tris base (121.14 g/mol)	242.5 g	2000 mM
EDTA (0.5 M)	200 mL of 0.5M EDTA	100 mM
Acetic acid (17.4 M)	57.45 mL	1000 mM
Adjust to pH 8.5 and dilute to 1x with Milli-Q water before use		