



Oct 16, 2018

Working

Manual dissection of the *Schistosoma mansoni* head and back end for transcriptomic analysis

Sarah Seyffert¹¹iGEM 2018 Heinrich-Heine-University[dx.doi.org/10.17504/protocols.io.uqnevve](https://doi.org/10.17504/protocols.io.uqnevve)

iGEM Duesseldorf 2018

 Sarah Seyffert 

ABSTRACT

This protocol describes how to prepare a 100x M2

1x M2 media could be created from the stock and supplemental stocks and trace metals could be added afterwards.

Recipes for standard and alternative M2 for culturing freshwater cyanobacteria, such as *Synechocystis* sp. PCC 6803, in a co-culture together with heterotrophic bacteria.

Media is usually not suitable for marine cyanobacteria.

Slightly modified from:

Weiss, Taylor L., Eric J. Young, and Daniel C. Ducat. "A synthetic, light-driven consortium of cyanobacteria and heterotrophic bacteria enables stable polyhydroxybutyrate production." *Metabolic engineering* 44 (2017): 236-245.

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

Important information

- 1 Always work under sterile conditions

Addition of ingredients for M2 in a 1 l bottle

- 2 Add $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ ($3.6 \text{ g} \cdot \text{L}^{-1}$)
 - If you do not have hydrated CaCl_2 , use 2.718 g CaCl_2 (powder) for 1 L of the stock.
 - 0.25 M $\text{Na}_2\text{-EDTA}$ stock: Dissolve 2.32 g $\text{Na}_2\text{-EDTA}$ (powder) and fill the bottle up with MilliQ (nuclease-free) water up to 25 mL. Afterwards, the pH was adjusted with diluted NaOH to pH = 8.0.
 -
 - K_2HPO_4 ($82,3 \text{ g} \cdot \text{L}^{-1}$)
- 3 Add Citric acid ($0.6 \text{ g} \cdot \text{L}^{-1}$)
- 4 Add NaNO_3 ($17 \text{ g} \cdot \text{L}^{-1}$)

- 5 Add $\text{MgSO}_4 \cdot 7 \text{H}_2\text{O}$ ($7.49 \text{ g} \cdot \text{L}^{-1}$)
 - If you do not have hydrated MgSO_4 , use 3.659 g MgSO_4 (powder) for 1 L of the stock.
- 6 $0.25 \text{ M Na}_2\text{-EDTA}$, pH 8.0 ($0.56 \text{ ml} \cdot \text{L}^{-1}$)
- 7 Fill up with MiliQ water to one liter
- 8 Mix the solution with a stir bar
- 9 autoclave afterwards



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited