

# Mouse P1 Kidney Cold-Active Protease Single Cell Dissociation

## Version 2

Andrew Potter, Steve Potter

### Abstract

Method used to derive single cell suspension from P1 mouse kidneys on ice, generating a cell suspension with greatly reduced artifact gene expression changes and suitable for downstream analysis using 10x Chromium or DropSeq.

**Citation:** Andrew Potter, Steve Potter Mouse P1 Kidney Cold-Active Protease Single Cell Dissociation. **protocols.io** dx.doi.org/10.17504/protocols.io.q7ddzi6

**Published:** 21 Jun 2018

## Guidelines

### Storage Conditions of Reagents

Reagent	Storage Condition
DPBS (ThermoFisher, 14190144)	4°C
0.5 M EDTA (Ambion, AM9260G)	room temp.
BSA (Sigma, A8806)	4°C
Protease from <i>Bacillus Licheniformis</i> (Sigma, P5380)	Store 100 µL aliquots (100 mg/mL) in DPBS at -80°C
DNase 1 (Applichem, A3778)	Store 10 µL aliquots (250 U/10 µL) in DPBS at -80°C

### Required Equipment

Equipment	Supplier	Catalog no.
gentleMACS dissociator	Miltenyi	130-093-235

### The protocol workflow is as follows:

- Isolate Kidney
- Initial digestion: triturate on ice
- Perform gentleMACS
- Continue triturating on ice
- Preparing cells for Chromium/DropSeq

BEFORE STARTING

**Prepare *Bacillus Licheniformis* enzyme mix just prior to starting dissociation:**

Volume (µl)	Reagent	Final concentration
894	DPBS	1X
1	0.5 M EDTA	0.5 mM
5	DNase 1 (250 U/10 µL)	125 U / mL

100

*B. Lich* (100 mg/mL)

10 mg/mL

**+25 mg tissue / 1 mL enzyme mix**

To prepare 0.01% BSA/PBS:

**Make stock of 10% BSA in DPBS and store at -20 °C. To make PBS/BSA 0.01% aliquot 50 mL of DPBS in 50 mL conical and pipet in 50 µL of 10% BSA stock.**

**Prepare 10% FBS/PBS with heat-inactivated FBS.**

## Protocol

### Step 1.

Extract & isolate P1 kidneys in ice-cold PBS.

---

### Step 2.

Mince kidneys on top of petri dish, on ice, using razor blade.

---

### Step 3.

Weigh out 25 mg of tissue for each tube of *B. Lich.* enzyme mix (2 tubes total).

 **AMOUNT**

25 mg Additional info:

---

### Step 4.

Incubate tissue + enzyme on ice for 7 minutes while triturating 15 strokes using 1 mL pipet every 2 minutes set to 700 µL - first with tip cut off.

---

### Step 5.

Monitor digestion by taking small aliquot and visualizing under scope (every 5 minutes).

---

### Step 6.

After 7 minutes, take the digest mix (combine the two tubes) and pipet into Miltenyi C-tube (placed on ice); take C-tube to gentleMACS placed in 4° cold room. Run program brain\_03 two times.

 **TEMPERATURE**

4 °C Additional info:

---

### Step 7.

After MACS, briefly quick spin the MACS tube (to 500 G) at 4 °C to ensure contents are in the bottom of the tube.

 **TEMPERATURE**

4 °C Additional info:

---

### Step 8.

Re-suspend and visualize cells using scope by taking small aliquot and using a slide; continue digesting cells in C-tube on ice for 8 additional minutes while triturating every 2 min 15 strokes using a 1 mL pipet.

---

### Step 9.

Add 3 mL ice-cold 10% FBS/PBS to digest mix in C-tube to inhibit the protease.

 [AMOUNT](#)

3 ml Additional info: ice-cold 10% FBS/PBS

---

### Step 10.

Transfer digest mix to a 15 mL conical. Spin 600 G for 5 minutes at 4 °C; discard supernatant; re-suspend cell pellet in 2 mL ice-cold PBS/BSA.

 [AMOUNT](#)

2 ml Additional info: re-suspend in PBS/BSA

 [TEMPERATURE](#)

4 °C Additional info:

---

### Step 11.

Filter re-suspended cells using 30 µm filter on sterile 50 mL conical on ice - rinse filter with 4 mL ice-cold PBS/BSA. Transfer flow-through to 15 mL conical.

 [AMOUNT](#)

4 ml Additional info: rinse filter

---

### Step 12.

Spin 15 mL conical tube containing filtered cells 600 G for 5 minutes at 4 °C; discard supernatant and re-suspend pellet in 10 mL ice-cold PBS/BSA.

 [AMOUNT](#)

10 ml Additional info: PBS/BSA

 [TEMPERATURE](#)

4 °C Additional info:

---

### Step 13.

Repeat rinse/spin in ice-cold PBS/BSA.

---

### Step 14.

Remove supernatant and re-suspend in 1-2 mL ice-cold PBS/BSA.

---

### Step 15.

Examine using hemocytometer and adjust concentration to 100 cells/uL for DropSeq or 1,000 cells/μL for 10X Chromium.

---