

# Counting Cells with Hemocytometer

Tyler C. Moore

## Abstract

Counting cells with a hemocytometer is an easy way to determine relatively accurate numbers of viable cells. After determining cell counts, cells can be passaged, frozen away, or used for an experiment at a particular density.

**Citation:** Tyler C. Moore Counting Cells with Hemocytometer. **protocols.io**

dx.doi.org/10.17504/protocols.io.nxsdfne

**Published:** 20 Mar 2018

## Guidelines

Keep cells sterile by wearing gloves, ethanol treating your hands and surfaces, using sterile pipette tips, and **only opening vials with cells in the BSL-2 laminar flow hood.**

## Materials

Trypan Blue 100 mL [7050](#) by [Stemcell Technologies](#)

Complete DMEM (DMEM 10% HI FBS 50ug/mL Gentamycin) by [Gibco - Thermo Fischer](#)

- ✓ Hemocytometer (Neubauer) by Contributed by users
- ✓ Compound Microscope by Contributed by users

## Protocol

### Step 1.

#### Collect cell suspension

Cell suspension must be in a known volume. Adherent cells can be detached by methods appropriate for your particular cell line.

### Step 2.

#### Dilute an aliquot of cells in Trypan Blue

Use a 1:10 dilution for many cells, 1:2 dilution if working with fewer cells



$\text{Cells/mL} = (\text{Cell count})/(\text{number of chambers counted}) \times \text{dilution} \times 10^4$

EX: If you diluted your cells 1:10 in trypan blue, counted chambers 1 and 4, and counted a total of 150 cells:  $150/2 = 75 \times 10 \times 10^4 = 7.5 \times 10^6$  cells/mL

NOTE: Concentration (cells/mL) x total volume (mL) = total number of cells