



Adult mouse ear dissociation (on ice)

Andrew Potter¹

¹CCHMC

dx.doi.org/10.17504/protocols.io.ucwesxe

Working

Oct 05, 2018

Human Cell Atlas Method Development Community



Andrew Potter



ABSTRACT

This protocol is used to dissociate adult (8-10 wk) mouse ear "on ice" (in order to reduce gene expression artifacts). The cell yield is 9000 cells/mg (213000 cells total from 23 mg tissue), with 98% viable. The protocol involves a 2-layered dissociation, incubating on ice with *Bacillus licheniformis* enzyme. The dissociation is augmented by the use of a dounce homogenizer to help break up the tissue.



Ears adult mouse
dissociation
protocol.pdf

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

***Bacillus Licheniformis* Enzyme Mix (1 mL per 23 mg tissue):**

100 μ L b. lich 100 mg/mL (10 mg/mL final conc. - Sigma, P5380)

1 μ L 0.5 M EDTA (0.5 mM final conc. - Sigma, A8806)

899 μ L DPBS (no Ca, Mg) ThermoFisher (cat. #14190)

Preparing enzymes:

The enzyme is made up in DPBS (#14190). It is aliquoted and stored at -80 °C at 100 mg/mL in 100 μ L aliquots.

Reagents

Bovine Serum Albumin - Sigma (A8806)

Hypothermosol FRS – Sigma (H4416)

Red Blood Cell Lysing Buffer – Sigma (R7757)

Required supplies:

2 mL dounce homogenizer – Bellco (1984-10002)

Centrifuge for 1.5 mL, 15 mL conicals

Pipettes and pipet tips

15 mL Conicals (MLS)

1.5 mL tubes (MLS)

30 μ M filters - Miltenyi (130-098-458)

Petri dishes (MLS)

Razor blades (MLS)

Ice bucket w/ice

Hemocytometers - InCyto Neubauer Improved (DHC-NO1-5)

Isolating tissue

- 1 Dissect out ear tissue and place in ice-cold hypothermosol on ice.
- 2 Using forceps, transfer ear tissue to petri dish on ice. **Mince ear tissue thoroughly on ice for 3-4 min on ice into 1-mm³ pieces** using razor blade while manipulating tissue with forceps – you will need to use grinding motion and vigorously break up tissue.

🕒 00:04:00 mince tissue

1st layer

- 3 Place **23 mg minced tissue** into 1 mL B. Lich enzyme cocktail. Incubate on ice. 🗂 23 mg tissue
- 4 Shake every min; triturate every 2 min with p1000 w/tip cut (start triturating at 2 min) for 20 minutes.
🕒 00:01:00 shake vigorously 🕒 00:02:00 triturate 10x 🕒 00:20:00 digest on ice
- 5 After 20 mins of triturating on ice, use pipet to **transfer digest mix to 2 mL dounce homogenizer**. Use 10 strokes of Pestle A every 2 min (4 series total, 8 min). Digest mix should become turbid. 🕒 00:02:00 dounce homogenize
- 6 At **28 min total digest time**, let tissue chunks settle in dounce homogenizer for 2 min. Transfer **90% of supernatant (containing released cells)** using 1 mL serological pipet onto a 30 μ M filter placed on a 15 mL conical.
🕒 00:02:00 let tissue settle
- 7 Rinse filter w/5 mL ice-cold PBS/BSA 0.04%. **Save flow through containing released cells on ice and keep filter on tube for 2nd layer.**

2nd layer

- 8 Add additional 1 mL b. lich enzyme mix to residual tissue chunks in dounce homogenizer. Continue incubating on ice.
- 9 **Use 10 strokes of pestle A every 2 min. for 15 additional minutes** (45 min. total time).
🕒 00:02:00 dounce 🕒 00:15:00 additional digest time
- 10 **Transfer entire volume to same 30 μ M filter on 15 mL conical**. Rinse with additional 5 mL ice-cold PBS/BA 0.04%.

Preparing cells for single cell sequencing

- 11 Centrifuge cells for 300 g for 5 min at 4 °C.
- 12 Remove supernatant & re-suspend in 100 μ L ice-cold PBS/BSA 0.04% (leave in 15 mL conical).
- 13 Add 1 mL RBC lysis buffer to cell mix. Triturate 10X. Let sit on ice for 2 min.
- 14 Add 12 mL ice-cold PBS/BSA 0.04% to dilute RBC lysis buffer.
- 15 Spin 300 g for five min. Remove supernatant. Re-suspend in 100-200 μ L PBS/BSA 0.04%.
- 16 Examine using hemocytometer with trypan blue.
- 17 Adjust concentration to 1,000 cells/ μ L for Chromium or 100 cells/ μ L for DropSeq.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited