



Feb 27, 2019

Working

Preparation of feeder-free iPSCs culture

In 1 collection

Celeste Karch¹, Rita Martinez¹, Jacob Marsh¹¹Washington University in St Louis

dx.doi.org/10.17504/protocols.io.x83fryn

Neurodegeneration Method Development Community

Tech. support email: ndcn-help@chanzuckerberg.com

Celeste Karch

Washington University in St Louis

IPSC CORTICAL
DIFFERENTIATION
022017.pdf

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

This protocol is part of the [IPSC CORTICAL DIFFERENTIATION](#) collection.

This method should be performed using sterile technique.

MATERIALS TEXT

Please refer to the attached full manuscript for required materials.

SAFETY WARNINGS

Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.

1 Coat appropriate sized tissue culture dish with Matrigel. Swirl plate to distribute evenly across surface of dish.












2 Allow Matrigel to set for 01:00:00 at room temperature or in a 37 °C humidified chamber.



This will provide enough cells for a full 96 well plate of neural aggregates.



We typically thaw 1 vial of iPSC into 1 well of a 6 well (1mL Matrigel per well). Volumes below are based on 6 well plates.

- 3 Prepare mTeSR1 feeder-free medium by aseptically adding  **100 ml** of thawed 5x supplement to  **400 ml** Basal Medium.
- 4 Add  **5 ml** penicillin/streptomycin if needed. Mix thoroughly and warm to room temperature.
- 5 Thaw iPSCs ($1-2 \times 10^6$ cells/mL) in a  **37 °C** water bath.
- 6 Add  **1 ml** of thawed iPSC to  **9 ml** of DMEM/F12 in a 15mL conical tube (Do not mix). Centrifuge iPSC at 750 rpm for  **00:03:00**.
- 7 Carefully aspirate medium from iPSC pellet. Add  **2 ml** of mTeSR1 supplemented with Rock inhibitor (10 μ M final) - **avoid breaking up clumps of suspended cells**.
- 8 Aspirate Matrigel from tissue culture dish. Add  **2 ml** of suspended iPSCs. Place plate in  **37 °C**, 5% CO₂ and 95% humidified chamber, making a “T” motion to ensure even distribution of cells.
- 9 Change medium with  **2 ml** of fresh mTesR1 medium daily.



As iPSC become more confluent, increase mTesR to 3mL per well.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited