

Jan 03,  
2020

## CODEX - Poly-L-Lysine Cover-Slip Preparation V.3

Leigh Propper<sup>1</sup>, Francesca Farris<sup>1</sup>, Marda Jorgensen<sup>1</sup>

<sup>1</sup>University of Florida

**1** Works for me dx.doi.org/10.17504/protocols.io.baxyifpw

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

### ABSTRACT

#### Poly-Lysine Coverslip Preparation

This protocol describes the process of creating Poly-lysine-coated coverslips that are required for mounting tissue sections for the CODEX® experiment workflow.

### GUIDELINES

- Managers and supervisors - are responsible for making sure that technicians are properly trained and equipment and facility are maintained in good working order.
- Laboratory personnel - are responsible for reading and understanding this SOP and related documents and to perform these tasks in accordance with the SOPs.

### MATERIALS

NAME	CATALOG #	VENDOR
WypAll™ L30 Multipurpose Wipers	19168203	Thermo Fisher
Parafilm Wrap PM996, 4 in. Wide, 125 ft./Roll	CP0672040	Thermo Fisher
Dumont Forceps (Cover Slip Forceps)	11251-33	Fine Science Tools
EMS Glass Cover Slips 22mm x 22mm	72204-01	Electron Microscopy Sciences
Poly-L-lysine 0.1% (w / v)	P8920	Sigma Aldrich

### SAFETY WARNINGS

- Use physical safety precautions when working with sharps (disposable blades, coverslips, etc).

### BEFORE STARTING

- Wear gloves when handling cover slips to avoid depositing epithelial cells or residues on surfaces.
- You will need a sanitized or acid-washed beaker for this procedure.

- 1 Remove 12-15 cover-slips from box. 1m
- Cover-Slip brand and type are REQUIRED, not suggested.**



CODEX requires use of these specific cover-slips for their diagnostic process.

- 2 Pour enough poly-l-lysine solution to cover the bottom of the clean beaker.
- 3 Gently place the desired amount of cover-slips into the beaker, spreading them out and allowing both surfaces to have contact with the solution. 1m
- 4 Slowly swirl the beaker to spread the cover-slips around the base. 30s
- 5 Add 7 mL of poly-lysine solution above the cover-slips to ensure that all are fully covered. 2m

#### 7 ml Poly-Lysine

- 6 Cover the beaker with plastic wrap or parafilm, and secure tightly with a rubber band to prevent evaporation. 2m
- 7 Leave cover-slips in poly-lysine solution for a **minimum** of 12 hours and up to one week **at room temperature**. 12h

 **12:00:00 or up to one week.**

- 8 After waiting the required incubation period, carefully remove the rubber band and parafilm from around the neck of the beaker. 2m
- 9 Gradually pour the remaining poly-lysine solution into the proper waste container. 2m
- 10 Fill the same beaker containing the cover-slips to half volume with double-distilled water(ddH<sub>2</sub>O) or purified water. 2m
- 11 Swirl the contents to mix the solution. 1m


- 12 Let the beaker and cover-slips sit for 30 seconds. 30s
- To prevent removal of poly-lysine, do not soak in water for >1 minute during each washing step.**
- 13 Slowly pour off the water into the sink. This completes **Wash #1** 30s
- 14 [go to step #10 and repeat steps 9-12 \( Wash #2\)](#) 1m 30s
- 15 [go to step #10 and repeat steps 9-12 \( Wash #3\)](#) 1m 30s
- 16 [go to step #10 and repeat steps 9-12 \( Wash #4\)](#) 1m 30s
- 17 [go to step #10 and repeat steps 9-12 \( Wash #5\)](#) 1m 30s
- 18 [go to step #10 and repeat steps 9-12 \( Wash #6\)](#) 1m 30s
- 19 [go to step #10 and repeat steps 9-12 \( Wash #7\)](#) 1m 30s
- 20 After completing the 7 washes, place 2 lint-free Wypall towels on the bench top. 30s
- 21 Remove the cover-slips from the water, placing them on top of the first set of towels. Ensure the cover-slips are not overlapping to allow proper drying. 2m
- 22 Using the specialty forceps indicated for use by CODEX, flip over each cover-slip onto the second clean towel to dry the reverse side. 20m



Specialty forceps indicated in the materials section.

- 23 Leave the cover-slips on the Wypall towel to dry. 24h

- 24 When the cover-slips are dry, the Poly-Lysine-coated cover-slips can be stored in a sterile petri dish or similarly covered container<sup>30m</sup> for up to 2 (two) months.



Sterile Petri Dish containing the finished Poly-L-Lysine coated cover-slips for CODEX processing.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited