

# A Typical DNase I Reaction (M0303)

## New England Biolabs

### Abstract


This is a protocol for a typical DNase I Reaction, using the M0303 RNase-free DNase I.

**Citation:** New England Biolabs A Typical DNase I Reaction (M0303). **protocols.io**

dx.doi.org/10.17504/protocols.io.cqevtd

**Published:** 27 Jan 2015

## Materials

 DNase I (RNase-free) - 1,000 units [M0303S](#) by [New England Biolabs](#)

## Protocol

### Step 1.

Resuspend **10 µg** RNA in 1X DNase I Reaction Buffer to a final volume of **100 µl**

### Step 2.

Add 2 units of DNase I

### Step 3.

Mix thoroughly

### Step 4.

Incubate at 37°C for 10 minutes.

 **DURATION**

00:10:00

### Step 5.

Add **1 µl** of 0.5 M EDTA (to a final concentration of **5 mM**).

 **AMOUNT**

1 µl Additional info:

 **DURATION**

00:10:00

 **ANNOTATIONS**

**Cecilie Christensen** 07 Mar 2018

Is the pH 8 for the EDTA?

### Step 6.

Heat inactivate at 75°C for 10 minutes.